




# QUALITY MEETS PERFORMANCE

## Kangdi KYQB Linear Cutter


# Kangdi KYQB Linear Cutter

## Product Features


- New knife blade with each firing
- Quick release button
- Safety mechanism when cartridge is fired



**Rear or Center Hinge Clamping**  
for one handed positioning




**Safety mechanism**  
to prevent firing  
spent reload



**Quick release button** for easy unclamping and removal of forks

**Disposable loading unit**  
with secure grip and new blade with each firing



Model	Reload	Color	Staple Height	Tissue Thickness
KYQB-60L	KYQZ-60L	Blue	3.8mm	1.5
KYQB-60H	KYQZ-60H	Green	4.5mm	2
KYQB-80L	KYQZ-80L	Blue	3.8mm	1.5
KYQB-80H	KYQZ-80H	Green	4.5mm	2
KYQB-100L	KYQZ-100L	Blue	3.5mm	1.5
KYQB-100H	KYQZ-100H	Green	4.5mm	2

For more information or to try a Kangdi product in your next procedure, email [customers@kangdi-med.com](mailto:customers@kangdi-med.com)

US150366