



SHAOXING
CARERE MEDICAL
APPLIANCE CO.,LTD



COMPANY PROFILE



Shaoxing CARERE Medical Appliance Co.,Ltd located in Pingxing Road, Pingshui Town, Keqiao District, Shaoxing City 312050, Zhejiang Province, China.

We are specialized in research, development, production and marketing of the Anesthesia and Respiratory Products. Our featured products are Anesthesia Breathing Circuit, Bacterial/Viral Filter (BV Filter), Heat Moisture Exchanger Filter(HMEF), Anesthetic Mask, Breathing Bag , Catheter Mount, Oxygen Mask, Laryngeal Mask Airway, Tracheal Tube and so on.The full varieties of products can meet different requests of customers. Our products are mainly exported to Europe, South America, South eastern Asia, Middle East, Africa and other countries and districts. Till now, we have set up the good long-term cooperation with our customers.

The company have 1000 square meters of 100,000-class clean room workshop equipped with more than 80 sets of advanced equipment for testing, manufacturing and packaging. We have established perfect quality management in accordance with the standards of National Medical Industry. Moreover, we have been ensured by a complete quality guarantee system of ISO 13485 with the products certified with CE.

With full confidence and absolute sincerity, the whole staff in CARERE company will adheres to the policy of "Quality and Efficiency, Customer First, Technology Innovation", to supply worldwide customers the best products and services.

CE ISO
0197 13485




CONTENTS

01	Anesthesia Circuits
04	Ventilator Circuits
07	Reusable Silicone Circuits
08	Catheter mount
11	Breathing System Filters
14	Anesthetic Masks
15	Oxygen Masks
18	Tracheal Tubes
20	Laryngeal Mask Airways
21	Breathing Bags
22	Tubing Series
23	Connectors&Accessories



ANESTHESIA CIRCUITS




Anesthesia Breathing Circuits are available corrugated, expandable and PVC smoothbore, which can be in combination with Anesthetic mask, Breathing system filter, Breathing bag and Respiratory machine, offer a simple, convenient and efficient passage for clinic gas delivery, such as anesthetic gas, oxygen gas. If you require a unique length or configuration, we can assemble a circuit to your exact specifications.


- 100% latex free
- Constructed of flexible and lightweight tubing
- High transparency for good visibility
- Wide variety of lengths accommodate all patient positioning needs
- Cuffed connections ensure good sealing
- Available in both adult and pediatric diameters

ANESTHESIA CIRCUIT–Corrugated


With 22M/15F gas sampling elbow, parallel-wye connector, corrugated tubes, end connectors, breathing bag



CODE	Configuration
CRR001AC	1.6m length, Adult
CRR002AC	1.5m length, Pediatric



CODE	Configuration
CRR003AC	1.6m length, Adult, 0.8m limb
CRR004AC	1.5m length, Pediatric, 0.75m limb




CODE	Configuration
CRR005AC	1.6m length, Adult, 0.8m limb, 2L breathing bag
CRR006AC	1.5m length, Pediatric, 0.75m limb, 1L breathing bag

*Special products can be manufactured upon request.



ANESTHESIA CIRCUITS

ANESTHESIA CIRCUIT–Smoothbore


With 22M/15F gas sampling elbow, parallel-wye connector, PVC smoothbore tubes, PVC end connectors, breathing bag



CODE	Configuration
CRR007AC	1.6m length, Adult
CRR008AC	1.5m length, Pediatric




CODE	Configuration
CRR009AC	1.6m length, Adult, 0.8m limb
CRR010AC	1.5m length, Pediatric, 0.75m limb




CODE	Configuration
CRR011AC	1.6m length, Adult, 0.8m limb, 2L breathing bag
CRR012AC	1.5m length, Pediatric, 0.75m limb, 1L breathing bag

ANESTHESIA CIRCUIT–Expandable


With 22M/15F gas sampling elbow, parallel-wye connector, expandable tubes, end connectors, breathing bag



CODE	Configuration
CRR013AC	1.6m length, Adult
CRR014AC	1.5m length, Pediatric



CODE	Configuration
CRR015AC	1.6m length, Adult, 0.8m limb
CRR016AC	1.5m length, Pediatric, 0.75m limb




CODE	Configuration
CRR017AC	1.6m length, Adult, 0.8m limb, 2L breathing bag
CRR018AC	1.5m length, Pediatric, 0.75m limb, 1L breathing bag

*Special products can be manufactured upon request.

ANESTHESIA CIRCUITS

ANESTHESIA CIRCUIT-Expandable Kit




CODE	Configuration
CRR019AC	1.8m length,Adult,3L breathing bag
CRR020AC	2.4m length, Adult, 3L breathing bag
CRR021AC	1.8m length, Adult, 3L breathing bag,filter
CRR022AC	2.4m length, Adult, 3L breathing bag,filter
CRR023AC	1.8m length,Adult,3L breathing bag,filter, size 5 anesthetic mask,gas sampling line



CODE	Configuration
CRR024AC	2.4m length,Adult,3L breathing bag,filter, size 5 anesthetic mask,gas sampling line
CRR025AC	1.8m length, Pediatric, 1L breathing bag,filter
CRR027AC	1.8m length,Pediatric,1L breathing bag,filter, size 3 anesthetic mask,gas sampling line

ANESTHESIA CIRCUIT-Coaxial Circuit




CODE	Configuration
CRR201BC	1.8m length,Adult
CRR202BC	1.8m length,Adult,2L breathing bag
CRR203BC	1.8m length,Adult,3L breathing bag
CRR204BC	1.8m length,Adult, 3L breathing bag,filter



CODE	Configuration
CRR205BC	1.8m length,Adult, 3L breathing bag,filter, size 5 anesthetic mask
CRR206BC	1.8m length,Adult, 3L breathing bag,filter, size 5 anesthetic mask,gas sampling line

ANESTHESIA CIRCUIT-Duo-limb Circuit



CODE	Configuration
CRR211DC	1.8m length,Adult
CRR212DC	1.8m length,Adult,2L breathing bag
CRR213DC	1.8m length,Adult,3L breathing bag
CRR214DC	1.8m length,Adult, 3L breathing bag,filter



CODE	Configuration
CRR215DC	1.8m length,Adult, 3L breathing bag,filter, size 5 anesthetic mask
CRR216DC	1.8m length,Adult, 3L breathing bag,filter, size 5 anesthetic mask,gas sampling line


*Special products can be manufactured upon request.

VENTILATOR CIRCUITS



VENTILATOR CIRCUIT-Corrugated


With 22M/15F gas sampling elbow, parallel-wye connector with ports, corrugated tubes, water traps, end connectors, hanger, gas sampling line, safety cap



CODE

Configuration


CRR101VC	1.6m length, Adult, 1 water trap
CRR102VC	1.5m length, Pediatric, 1 water trap



CODE

Configuration


CRR103VC	1.6m length, Adult, 1 water trap, 0.8m limb
CRR104VC	1.5m length, Pediatric, 1 water trap, 0.75m limb



CODE

Configuration


CRR105VC	1.6m length, Adult, 2 water traps
CRR106VC	1.5m length, Pediatric, 2 water traps



CODE

Configuration


CRR107VC	1.6m length, Adult, 2 water traps, 0.8m limb
CRR108VC	1.5m length, Pediatric, 2 water traps, 0.75m limb



CODE

Configuration

CRR121VC	1.6m length, Adult, 2 water traps, 0.8m limb, hanger, gas sampling line
CRR122VC	1.5m length, Pediatric, 2 water traps, 0.75m limb, hanger, gas sampling



CODE

Configuration

CRR123VC	1.6m length, Adult, 2 water traps, 0.8m limb, hanger, safety cap
CRR124VC	1.5m length, Pediatric, 2 water traps, 0.75m limb, hanger, safety cap


*Special products can be manufactured upon request.

VENTILATOR CIRCUITS




VENTILATOR CIRCUIT-Smoothbore


With 22M/15F gas sampling elbow, parallel-wye connector with ports, PVC smoothbore tubes, water traps, PVC end connectors, hanger, safety cap




CODE	Configuration
CRR109VC	1.6m length, Adult, 1 water trap
CRR110VC	1.5m length, Pediatric, 1 water trap




CODE	Configuration
CRR111VC	1.6m length, Adult, 1 water trap, 0.8m limb
CRR112VC	1.5m length, Pediatric, 1 water trap, 0.75m limb




CODE	Configuration
CRR113VC	1.6m length, Adult, 2 water traps
CRR114VC	1.5m length, Pediatric, 2 water traps



CODE	Configuration
CRR115VC	1.6m length, Adult, 2 water traps, 0.8m limb
CRR116VC	1.5m length, Pediatric, 2 water traps, 0.75m limb



CODE	Configuration
CRR125VC	1.6m length, Adult, 2 water traps, 0.8m limb, safety cap
CRR126VC	1.5m length, Pediatric, 2 water traps, 0.75m limb, safety cap



CODE	Configuration
CRR127VC	1.6m length, Adult, 2 water traps, 0.8m limb, hanger, safety cap
CRR128VC	1.5m length, Pediatric, 2 water traps, 0.75m limb, hanger, safety cap


*Special products can be manufactured upon request.

VENTILATOR CIRCUITS



VENTILATOR CIRCUIT-Exhalation valve

With exhalation valve and line, nebulizer, pressure monitoring line, 22M/15F gas sampling elbow, corrugated/PVC smoothbore tubes, water trap, catheter mount, safety cap




CODE

Configuration

CRR131VC

1.6m length, exhalation valve and line,
15cm corrugated catheter mount




CODE

Configuration

CRR132VC

1.6m length, exhalation valve and line,
15cm smoothbore catheter mount




CODE

Configuration

CRR133VC

1.6m length, exhalation valve and line,
nebulizer, 15cm corrugated catheter mount




CODE

Configuration

CRR134VC

1.6m length, exhalation valve and line,
nebulizer, 15cm smoothbore catheter mount




CODE

Configuration

CRR135VC

1.8m length, Adult, ventilator circuit with
external exhalation valve



CODE



Configuration

CRR136VC


1.8m length, Adult, ventilator circuit with
external exhalation valve, water trap

*Special products can be manufactured upon request.


REUSABLE SILICONE CIRCUITS




CODE	Configuration
CRR051AC	1.5m length, Adult
CRR052AC	1.5m length, Pediatric
CRR053AC	1.5m length, Infant



CODE	Configuration
CRR054AC	1.5m length, Adult, 0.6m limb
CRR055AC	1.5m length, Pediatric, 0.6m limb
CRR056AC	1.5m length, Infant, 0.6m limb



CODE	Configuration
CRR141VC	1.5m length, Adult, 1 water trap, 0.6m limb
CRR142VC	1.5m length, Pediatric, 1 water trap, 0.6m limb
CRR143VC	1.5m length, Infant, 1 water trap, 0.6m limb




CODE	Configuration
CRR144VC	1.5m length, Adult, 2 water traps, 0.6m limb
CRR145VC	1.5m length, Pediatric, 2 water traps, 0.6m limb
CRR146VC	1.5m length, Infant, 2 water traps, 0.6m limb

*Special products can be manufactured upon request.

CATHETER MOUNT

CATHETER MOUNT-Double Swivel Elbow

Available with corrugated, expandable and PVC smoothbore tubes. Transparency patient connectors provide visibility during surgery. A variety of connectors provide more options for airway management. ISO connections for a secure fit.




CODE

CRR301CM

Configuration

Φ15mm, 15cm corrugated tubing with double swivel elbow, 22F connector




CODE

CRR302CM

Configuration

Φ15mm, 15cm corrugated tubing with double swivel elbow, 15M connector




CODE

CRR303CM

Configuration

Φ15mm, 15cm smoothbore tubing with double swivel elbow, 22F connector




CODE

CRR304CM

Configuration

Φ15mm, 15cm smoothbore tubing with double swivel elbow, 15M connector




CODE

CRR305CM

Configuration

Φ15mm, 15cm expandable tubing with double swivel elbow, 22F connector



CODE

CRR306CM


Configuration

Φ15mm, 15cm expandable tubing with double swivel elbow, 15M connector

*Special products can be manufactured upon request.

CATHETER MOUNT

CATHETER MOUNT-Gas Sampling Elbow




CODE

Configuration

CRR311CM

Φ15mm,15cm corrugated tubing with 22M/15F gas sampling elbow,22F connector




CODE

Configuration

CRR312CM

Φ15mm,15cm corrugated tubing with 22M/15F gas sampling elbow ,15M connector




CODE

Configuration

CRR313CM

Φ15mm,15cm smoothbore tubing with 22M/15F gas sampling elbow,22F connector




CODE

Configuration

CRR314CM

Φ15mm,15cm smoothbore tubing with 22M/15F gas sampling elbow,15M connector




CODE

Configuration

CRR315CM

Φ15mm,15cm expandable tubing with 22M/15F gas sampling elbow,22F connector



CODE


Configuration

CRR316CM


Φ15mm,15cm expandable tubing with 22M/15F gas sampling elbow,15M connector

CATHETER MOUNT


CATHETER MOUNT-Plain Elbow




CODE	Configuration
CRR351CM	Φ15mm,15cm corrugated tubing with 22M/15F plainelbow,22F connector




CODE	Configuration
CRR352CM	Φ15mm,15cm corrugated tubing with 22M/15F plainelbow , 15M connector




CODE	Configuration
CRR353CM	Φ15mm,15cm smoothbore tubing with 22M/15F plainelbow,22F connector



CODE	Configuration
CRR354CM	Φ15mm,15cm smoothbore tubing with 22M/15F plainelbow , 15M connector



CODE	Configuration
CRR355CM	Φ15mm,15cm expandable tubing with 22M/15F plainelbow,22F connector



CODE	Configuration
CRR356CM	Φ15mm,15cm expandable tubing with 22M/15F plainelbow , 15M connector

BREATHING SYSTEM FILTERS

BREATHING SYSTEM FILTERS



We provide an extensive line of Breathing System filter that meet all ISO standards and offer a high level of clinical protection to patients of any size.

Hygroscopic membrane traps patient's moisture providing effective airway humidification.

- Designed in light weight, enhancing patient's comfort
- Gas sampling ports available with tethered cap or luer cap
- Pediatric and adult size available

Bacteria Filters



CODE	CRR401BF	CRR402BF	CRR403BF
Description	BVF-01	BVF-02	BVF-03
Dead space	30ml	30ml	35ml
Pressure drop@30LPM	60Pa	60Pa	60Pa
Moisture output	/	/	/
Bacterial filtration efficiency	99.999%	99.999%	99.999%
Viral filtration efficiency	99.999%	99.999%	99.999%
Weight	23g	22g	22g
connections	22M/15F-22F/15M	22M/15F-22F/15M	22M/15F-22F/15M
Gas monitoring	Yes	No	Yes

HEPA



CODE	CRR401HP
Description	HEPA-01
Dead space	60ml
Tidal volume range	150-1500ml
Pressure drop@30LPM	30Pa
Bacterial filtration efficiency	99.9999%
Viral filtration efficiency	99.9999%
Weight	46g
connections	22M/15F-22F/15M
Gas monitoring	Yes

BREATHING SYSTEM FILTERS

HME Trach



CODE	CRR401HT	CRR402HT	CRR403HT	CRR404HT
Description	HME Trach basic with foam	HME Trach basic with paper	HME Trach double with foam	HME Trach double with paper
Dead space	16ml	16ml	16ml	16ml
Pressure drop@30LPM	50Pa	50Pa	50Pa	50Pa
Moisture output	24mg@500VT	24mg@500VT	24mg@500VT	24mg@500VT
Tidal volume range	200-1000ml	200-1000ml	200-1000ml	200-1000ml
weight	5g	9g	7g	5g
connections	15F	15F	15F	15F

Spirometry Filter



CODE	CRR401SF
Description	Spirometry filter
Dead space	75ml
Media type	Electrostatic
Bacterial filtration efficiency	99.999%
Viral filtration efficiency	99.999%
Pressure drop@30LPM	20Pa
Option	Mouth piece
Weight	48g
Connections	30F/25M

HMEF Series



CODE	CRR401HF	CRR402HF	CRR405HF	CRR406HF	CRR407HF
Description	HMEF Pediatric	HMEF Adult	HMEF Adult	HMEF Adult	HMEF Adult large
Material Type	K-resin + paper	K-resin + paper	PP +paper	PP +paper	PP +paper
Dead space	12ml	40ml	35ml	40ml	60ml
Tidal volume range	150-300ml	150-1500ml	150-1500ml	150-1500ml	300-1500ml
Pressure drop@30LPM	120Pa	100Pa	100Pa	100Pa	100Pa
Moisture output	24mg@500VT	32mg@500VT	32mg@500VT	32mg@500VT	34mg@500VT
Bacterial filtration efficiency	99.999%	99.999%	99.999%	99.999%	99.999%
Viral filtration efficiency	99.999%	99.999%	99.999%	99.999%	99.999%
Weight	15g	30g	25g	25g	31g
connections	22M/15F-22F/15M	22M/15F-22F/15M	22M/15F-22F/15M	22M/15F-22F/15M	22M/15F-22F/15M
Gas monitoring	Yes	Yes	Yes	Yes	Yes

BREATHING SYSTEM FILTERS

HMEF Series



CODE	CRR403HF	CRR404HF	CRR408HF
Description	HMEF Pediatric	HMEF Adult	HMEF Adult
Material Type	K-resin + foam	K-resin + foam	PP +foam
Dead space	12ml	40ml	35ml
Tidal volume range	150-300ml	150-1500ml	150-1500ml
Pressure drop@30LPM	120Pa	100Pa	100Pa
Moisture output	24mg@500VT	32mg@500VT	32mg@500VT
Bacterial filtration efficiency	99.999%	99.999%	99.999%
Viral filtration efficiency	99.999%	99.999%	99.999%
Weight	15g	25g	25g
connections	22M/15F-22F/15M	22M/15F-22F/15M	22M/15F-22F/15M
Gas monitoring	Yes	Yes	Yes

Salted HMEF Series



CODE	CRR409HF	CRR410HF
Description	HMEF Pediatric	HMEF Adult
Material Type	K-resin + salted paper	K-resin + salted paper
Dead space	12ml	40ml
Tidal volume range	150-300ml	150-1500ml
Pressure drop@30LPM	120Pa	100Pa
Moisture output	30mg@500VT	38mg@500VT
Bacterial filtration efficiency	99.999%	99.999%
Viral filtration efficiency	99.999%	99.999%
Weight	15g	31g
connections	22M/15F-22F/15M	22M/15F-22F/15M
Gas monitoring	Yes	Yes

HMEF with Catheter Mount



CODE	Configuration
CRR421HC	Adult HMEF+15cm expandable tubing with 22F connector & 22M/15F connector
CRR422HC	Adult HMEF+15cm expandable tubing with 15M connector & 22M/15F connector
CRR423HC	Pediatric HMEF+15cm expandable tubing with 22F connector & 22M/15F connector
CRR424HC	Pediatric HMEF+15cm expandable tubing with 15M connector & 22M/15F connector



CODE	Configuration
CRR425HC	Adult HMEF + 15cm expandable tubing with 22F connector & 22M/15F connector, red cap
CRR426HC	Adult HMEF + 15cm expandable tubing with 15M connector & 22M/15F connector, red cap
CRR427HC	Pediatric HMEF + 15cm expandable tubing with 22F connector & 22M/15F connector, red cap
CRR428HC	Pediatric HMEF + 15cm expandable tubing with 15M connector & 22M/15F connector, red cap

ANESTHETIC MASKS

Anesthetic Masks

- 100% latex-free
- ISO standard connectors for secure connection
- Inflation valve for easy cushion pressure adjustment
- Soft air cushion for good sealing and patient comfort
- Transparent cone for excellent visibility

Anesthetic Mask with Valve Horizontal



CODE	Configuration
CRR501AM	Size-0, Neonate,
CRR502AM	Size-1, Infant,
CRR503AM	Size-2, Pediatric, yellow hook
CRR504AM	Size-3, Adult small, green hook
CRR505AM	Size-4, Adult medium, red hook
CRR506AM	Size-5, Adult large, blue hook

Anesthetic Mask with Valve Up



CODE	Configuration
CRR507AM	Size-0, Neonate
CRR508AM	Size-1, Infant
CRR509AM	Size-2, Pediatric, yellow hook
CRR510AM	Size-3, Adult small, pink hook
CRR511AM	Size-4, Adult medium, green hook
CRR512AM	Size-5, Adult large, blue hook

Anesthetic Mask with Non-valve



CODE	Configuration
CRR513AM	Size-0, Neonate, blue
CRR514AM	Size-1, Infant, orange
CRR515AM	Size-2, Pediatric, pink
CRR516AM	Size-3, Adult small, white
CRR517AM	Size-4, Adult medium, yellow
CRR518AM	Size-5, Adult large, green

Silicone Anesthetic Mask

- latex-free and reusable
- Transparent for visual check of bleeding, vomiting and spontaneous breathing
- All models are reusable and can be autoclaved repeatedly at 134°C



CODE	Configuration
CRR521AM	Size-0, Neonate
CRR522AM	Size-1, Infant
CRR523AM	Size-2, Pediatric
CRR524AM	Size-3, Adult small
CRR525AM	Size-4, Adult medium
CRR526AM	Size-5, Adult large
CRR527AM	Size-6, Adult extra large

OXYGEN MASKS



CODE	Configuration
CRR601OM	XL(Adult Elongated)
CRR602OM	L(Adult Asian)
CRR603OM	M(Pediatric Elongated)
CRR604OM	S(Pediatric Standard)

Oxygen Mask

- Made of clear, medical grade PVC
- With adjustable nose clip, elastic strap and latex free
- Medium concentration, approx 50%
- Color: green transparent, white transparent and light blue transparent is available.
- Length: with 2m tubing, different length can be customized
- Anit-crush tube can ensure oxygen follow even if the tube is kinked



CODE	Configuration
CRR605OM	XL(Adult Elongated)
CRR606OM	L(Adult Asian)
CRR607OM	M(Pediatric Elongated)
CRR608OM	S(Pediatric Standard)

Oxygen Mask with Reservoir Bag

- Oxygen Mask with Reservoir Bag
- Made of clear, medical grade PVC
- With reservoir bag
- With adjustable nose clip, elastic strap and latex free
- Color: green transparent, white transparent and light blue transparent is available
- Length: with 2m tubing, different length can be customized
- Anit-crush tube can ensure oxygen follow even if the tube is kinked

OXYGEN MASKS



CODE	Configuration
CRR609OM	XL(Adult Elongated)
CRR610OM	L(Adult Asian)
CRR611OM	M(Pediatric Elongated)
CRR612OM	S(Pediatric Standard)

Adjusted Oxygen Concentration	Oxygen flow
24%	3L/min
26%	3L/min
28%	6L/min
30%	6L/min
35%	9L/min
40%	12L/min
50%	15L/min

Venturi Mask

- Made of clear, medical grade PVC
- With adjustable nose clip, elastic strap and latex free
- Two color coded diluters for seven concentrations
- Green diluter for 24%,26%,28%,30%,white diluter for 35%,40%,50%
- Color: green transparent, white transparent and light blue transparent is available.
- Length: with 2m tubing, different length can be customized
- Anit-crush tube can ensure oxygen follow even if the tube is kinked



CODE	Configuration
CRR613OM	XL(Adult Elongated)
CRR614OM	L(Adult Asian)
CRR615OM	M(Pediatric Elongated)
CRR616OM	S(Pediatric Standard)

Adjusted Oxygen Concentration	Oxygen flow
24%(blue)	4L/min
28%(yellow)	4L/min
31%(white)	6L/min
35%(dark green)	8L/min
40%(pink)	8L/min
50%(orange)	10L/min

Adjustable Venturi Mask

- Made of clear, medical grade PVC
- With adjustable nose clip, elastic strap and latex free
- Single diluter for one concentration, six diluters can be selected
- Color: green transparent, white transparent and light blue transparent is available
- Length: with 2m tubing, different length can be customized
- Anit-crush tube can ensure oxygen follow even if the tube is kinked

OXYGEN MASKS



CODE	Configuration
CRR617OM	XL(Adult Elongated)
CRR618OM	L(Adult Asian)
CRR619OM	M(Pediatric Elongated)
CRR620OM	S(Pediatric Standard)

Nebulizer Mask

- Made of clear, medical grade PVC
- Operates even at horizontal position capacity:6cc
- Disinfected with boiling water (only for jar)
- Different type of nebulizer jars are available
- With adjustable nose clip, elastic strap and latex free
- Color: green transparent, white transparent and light blue transparent is available
- Length: with 2m tubing, different length can be customized
- Anit-crush tube can ensure oxygen follow even if the tube is kinked



CODE	Configuration
CRR621OM	Curved Prong,S(Pediatric)
CRR622OM	Curved Prong,L(Adult)
CRR623OM	Straight Prong,XXS(Neonate)
CRR624OM	Straight Prong,XS(Infant)
CRR625OM	Straight Prong,S(Pediatric)
CRR626OM	Straight Prong,L(Adult)
CRR627OM	Soft Prong,XXS(Neonate)
CRR628OM	Soft Prong,XS(Infant)
CRR629OM	Soft Prong,S(Pediatric)
CRR630OM	Soft Prong,L(Adult)

Nasal Oxygen Cannula

- Made of clear, medical grade PVC
- Soft prong, Straight prong, curved prong is available
- Color: green transparent, white transparent and light blue transparent is available.
- Length: with 2m tubing, different length can be customized
- Anit-crush tube can ensure oxygen follow even if the tube is kinked


TRACHEAL TUBES



Endotracheal Tube(Oral/Nasal)

- Made of thermosensitive, clear, medical grade PVC
- Available both with cuff and without cuff
- With high-volume, low-pressure cuff
- Murphy eye to avoid complete respiratory obstruction
- Radiopaque line throughout the tube for X-ray visualization

CODE	SIZE With Cuff	CODE	SIZE W/O Cuff
CRR701ET	3.0mm	CRR751ET	2.5mm
CRR702ET	3.5mm	CRR752ET	3.0mm
CRR703ET	4.0mm	CRR753ET	3.5mm
CRR704ET	4.5mm	CRR754ET	4.0mm
CRR705ET	5.0mm	CRR755ET	4.5mm
CRR706ET	5.5mm	CRR756ET	5.0mm
CRR707ET	6.0mm	CRR757ET	5.5mm
CRR708ET	6.5mm	CRR758ET	6.0mm
CRR709ET	7.0mm	CRR759ET	6.5mm
CRR710ET	7.5mm	CRR760ET	7.0mm
CRR711ET	8.0mm	CRR761ET	7.5mm
CRR712ET	8.5mm	CRR762ET	8.0mm
CRR713ET	9.0mm	CRR763ET	8.5mm
CRR714ET	9.5mm	CRR764ET	9.0mm
CRR715ET	10.0mm	CRR765ET	9.5mm




Reinforced Endotracheal Tube(Oral/Nasal)

- Made of thermosensitive, clear, medical grade PVC
- Available both with cuff and without cuff
- With high-volume, low-pressure cuff
- Murphy eye to avoid complete respiratory obstruction
- Radiopaque line throughout the tube for X-ray visualization
- A stainless steel spring is inserted in the tube to minimize the risk of kinking or crushing.

CODE	SIZE With Cuff	CODE	SIZE W/O Cuff
CRR716ET	3.0mm	CRR766ET	2.5mm
CRR717ET	3.5mm	CRR767ET	3.0mm
CRR718ET	4.0mm	CRR768ET	3.5mm
CRR719ET	4.5mm	CRR769ET	4.0mm
CRR720ET	5.0mm	CRR770ET	4.5mm
CRR721ET	5.5mm	CRR771ET	5.0mm
CRR722ET	6.0mm	CRR772ET	5.5mm
CRR723ET	6.5mm	CRR773ET	6.0mm
CRR724ET	7.0mm	CRR774ET	6.5mm
CRR725ET	7.5mm	CRR775ET	7.0mm
CRR726ET	8.0mm	CRR776ET	7.5mm
CRR727ET	8.5mm	CRR777ET	8.0mm
CRR728ET	9.0mm	CRR778ET	8.5mm
CRR729ET	9.5mm	CRR779ET	9.0mm
CRR730ET	10.0mm	CRR780ET	9.5mm

TRACHEAL TUBES



CODE	SIZE With Cuff	CODE	SIZE W/O Cuff
CRR731ET	3.0mm	CRR781ET	3.0mm
CRR732ET	3.5mm	CRR782ET	3.5mm
CRR733ET	4.0mm	CRR783ET	4.0mm
CRR734ET	4.5mm	CRR784ET	4.5mm
CRR735ET	5.0mm	CRR785ET	5.0mm
CRR736ET	5.5mm	CRR786ET	5.5mm
CRR737ET	6.0mm	CRR787ET	6.0mm
CRR738ET	6.5mm	CRR788ET	6.5mm
CRR739ET	7.0mm	CRR789ET	7.0mm
CRR740ET	7.5mm	CRR790ET	7.5mm
CRR741ET	8.0mm	CRR791ET	8.0mm
CRR742ET	8.5mm	CRR792ET	8.5mm
CRR743ET	9.0mm	CRR793ET	9.0mm

Silicone Endotracheal Tube(Oral/Nasal)

- Made of medical grade silicone.
- Available both with cuff and without cuff
- Endotracheal Tube features softer and less evasive tip to prevent damage to tracheal tissue.
- With high-volume,low-pressure cuff
- A stainless steel spring is inserted in the tube to minimize the risk of kinking or crushing.

Double Lumen Endobronchial Tube



CODE	Configuration
CRR744ET	Fr 28 (250-28L,250-28R)
CRR745ET	Fr32(250-32L,250-32R)
CRR746ET	Fr35(250-35L,250-35R)
CRR747ET	Fr37(250-37L,250-37R)
CRR748ET	Fr39(250-39L,250-39R)
CRR749ET	Fr41(250-41L,250-41R)

Intubation Stylet



CODE	Configuration
CRR795ET	6# (2.5-4.5mm)
CRR796ET	10#(4.0-6.0mm)
CRR797ET	12#(5.5-6.5mm)
CRR798ET	14#(6.0-10.0mm)

LARYNGEAL MASK AIRWAYS



Silicone LMA

- Made of medical grade silicone
- Silicone cuff is made by one time, provides good and soft sealing
- Available also for reinforced flexible version
- Cuff color can be customized: flesh, blue and transparent colors available
- MRI compatible type is available

CODE	Configuration Silicone
CRR801LM	1.0(0-5kgs,4ml)
CRR802LM	1.5(5-10kgs,7ml)
CRR803LM	2.0(10-20kgs,10ml)
CRR804LM	2.5(20-30kgs,14ml)
CRR805LM	3.0(30-50kgs,20ml)
CRR806LM	4.0(50-70kgs,30ml)
CRR807LM	5.0(70-100kgs,40ml)

CODE	Configuration Reinforced Silicone
CRR831LM	1.0(0-5kgs,4ml)
CRR832LM	1.5(5-10kgs,7ml)
CRR833LM	2.0(10-20kgs,10ml)
CRR834LM	2.5(20-30kgs,14ml)
CRR835LM	3.0(30-50kgs,20ml)
CRR836LM	4.0(50-70kgs,30ml)
CRR837LM	5.0(70-100kgs,40ml)



PVC LMA

- Made of medical grade PVC
- Available also for reinforced flexible version
- Smooth cuff provides good and soft sealing
- Color coded, easy to identify sizes

CODE	Configuration PVC
CRR808LM	1.0(0-5kgs,4ml)
CRR809LM	1.5(5-10kgs,7ml)
CRR810LM	2.0(10-20kgs,10ml)
CRR811LM	2.5(20-30kgs,14ml)
CRR812LM	3.0(30-50kgs,20ml)
CRR813LM	4.0(50-70kgs,30ml)
CRR814LM	5.0(70-100kgs,40ml)

CODE	Configuration Reinforced PVC
CRR871LM	1.0(0-5kgs,4ml)
CRR872LM	1.5(5-10kgs,7ml)
CRR873LM	2.0(10-20kgs,10ml)
CRR874LM	2.5(20-30kgs,14ml)
CRR875LM	3.0(30-50kgs,20ml)
CRR876LM	4.0(50-70kgs,30ml)
CRR877LM	5.0(70-100kgs,40ml)

BREATHING BAGS

Breathing Bag



CODE	Configuration
CRR901BB	0.5 Liter, latex, blue
CRR902BB	1.0 Liter, latex, blue
CRR903BB	2.0 Liters, latex, blue
CRR904BB	3.0 Liters, latex, blue



CODE	Configuration
CRR905BB	0.5 Liter, latex-free, green
CRR906BB	1.0 Liter, latex-free, green
CRR907BB	2.0 Liters, latex-free, green
CRR908BB	3.0 Liters, latex-free, green



CODE	Configuration
CRR909BB	0.5 Liter, Silicone, white
CRR910BB	1.0 Liter, Silicone, white
CRR911BB	2.0 Liters, Silicone, white
CRR912BB	3.0 Liters, Silicone, white

TUBING SERIES

Segmented Coils



CODE	Configuration
CRR9011P	Φ22mm, 100'(30m), clear, 6"(15cm) segment length
CRR9012P	Φ22mm, 100'(30m), blue, 6"(15cm) segment length
CRR9013P	Φ15mm, 200'(60m), clear, 6"(15cm) segment length

CPAP Tubings




CODE	Configuration
CRR9024P	Φ22mm, 1.8m length, 22mm ID cuff
CRR9022P	Φ22mm, 1.2m length, 22mm ID cuff

Silicone Tubings



CODE	Configuration
CRR9031P	Φ22mm, 1.5m length, 22mm ID cuff
CRR9032P	Φ22mm, 0.6m length, 22mm ID cuff
CRR9033P	Φ15mm, 1.5m length, 22mm ID cuff
CRR9034P	Φ15mm, 0.6m length, 22mm ID cuff
CRR9035P	Φ10mm, 1.5m length, 22mm ID cuff
CRR9036P	Φ10mm, 0.6m length, 22mm ID cuff

PVC Smoothbore Tubings



CODE	Configuration
CRR9071P	Φ22mm, 0.8m length with cuff
CRR9072P	Φ22mm, 1.0m length with cuff
CRR9073P	Φ22mm, 1.6m length with cuff
CRR9075P	Φ15mm, 0.8m length with cuff
CRR9076P	Φ15mm, 1.2m length with cuff
CRR9077P	Φ15mm, 1.5m length with cuff

CONNECTORS&ACCESSORIES



Elbow Connectors



CRR8001P



CRR8002P



CRR8003P



CRR8005P



CRR8004P

Caps



CRR8011P



CRR8012P



CRR8013P



CRR8014P



CRR8015P

CONNECTORS&ACCESSORIES

End Connectors



CRR8021P



CRR8022P



CRR8023P



CRR8024P



CRR8025P



CRR8026P



CRR8027P

Wye Connectors



CRR8031P



CRR8032P



CRR8033P



CRR8034P



CRR8035P



CRR8036P



CRR8037P

Water Trap



CRR7001P

Water Trap



CRR7002P

Exhalation Valve



CRR7003P

Exhalation Valve



CRR7004P

Nebulizer



CRR7005P

Tube Hanger



CRR7006P

Tube Clip



CRR7007P

T-piece Connector Gas Sampling Line



CRR7008P



CRR7009P



DO OUR *BEST* FOR *HEALTH*

SHAOXING CARERE MEDICAL APPLIANCE CO.,LTD

Add: Pingxing Road,Loujia Village,Pingshui Town, Keqiao District,
Shaoxing City 312050, Zhejiang Province, China.

Tel: +86 - 575 - 85020569

Fax: +86 - 575 - 85090579

Email: Info@crrmedical.com

[Http:// www.crrmedical.com](http://www.crrmedical.com)

绍兴柯锐医疗用品股份有限公司

地址：浙江省绍兴市柯桥区平水镇娄家村平兴路

邮编：312050

电话：0575 - 85020569

传真：0575 - 85090579

邮箱：Info@crrmedical.com

网址：www.crrmedical.com

