

# Contents

**01-15**

## **DIAGNOSTIC**

Stethoscope  
Stethoscope & Neurological hammer  
Sphygmomanometer  
ECG products  
ENT diagnostic set  
Thermometer

**16-17**

## **GYNAECOLOGY**

**18-23**

## **DRESSING**

Plaster  
Cotton products  
Bandage  
Non-Woven products

**24-38**

## **HEALTH-CARE**

Surgical products  
Tube & Bag  
Gloves  
Injection products  
Rubber & Plastic  
Scales  
Hospital bed

**39-44**

## **HEALING**

Wheelchair  
Walkers  
Stick

**45-57**

## **EQUIPMENT**

Operating table  
Operating lamp  
X-Ray unit  
Dental products  
Sterilizer  
Instrument  
Suction unit

**58-62**

## **LABORATORY**

Microscope  
Laboratory products  
Others

# DIAGNOSTIC





### DT010 SINGLE HEAD STETHOSCOPE

Anodized aluminum chestpiece of diameter 46mm for adult, brass chrome plated binaural, white plastic or soft materials ear tips, and one Y-PVC tubing in assorted colors: Black, Yellow, Red, Pink, Blue, Orange, Green, White, Gray, Navy, Brown, Peach, Lavender, etc.



### DT011 SINGLE HEAD STETHOSCOPE(PEDIATRIC)

Anodized aluminum chestpiece of diameter 38mm for pediatric, brass chrome plated binaural, white plastic or soft materials ear tips, and one Y-PVC tubing in as sorted colors.



### DT012 SINGLE HEAD STETHOSCOPE WITH NON-CHILL RING

Anodized aluminum chestpiece matched non-chill rings with a wide range of tubing colors for choice. Diameter of head: 46mm, for adult use.



### DT013 SINGLE HEAD STETHOSCOPE WITH NON-CHILL RING(PEDIATRIC)

Anodized aluminum chestpiece matched non-chill rings with a wide range of tubing colors for choice. Diameter of head: 38mm, for pediatric use.



### DT014 SHARP NEEDLE TYPE STETHOSCOPE

Brass anodized head, standard tubing, diameter of head: 48mm



### DT015 SHARP NEEDLE TYPE STETHOSCOPE

Brass anodized head, Black chrome coated, standard tubing, diameter of head: 48mm

# STETHOSCOPE



### DT110 DUAL HEAD STETHOSCOPE

Anodized aluminum dual chestpiece of diameter 46mm for adult, brass chrome plated binaural, white plastic or soft materials ear tips. Black flexible and long-lasting PVC tubing. Attractive colored chestpiece matching colorful tubing for your choice



### DT111 DUAL HEAD STETHOSCOPE(PEDIATRIC)

Anodized aluminum dual chestpiece of diameter 38mm for pediatric, brass chrome plated binaural, white plastic or soft materials ear tips, Black flexible and long-lasting PVC tubing. Attractive colored chestpiece matching color tubing for your choice.



### DT112 DUAL HEAD STETHOSCOPE WITH NON-CHILL RING

Anodized aluminum dual chestpiece matched non-chill rings with a wide range of tubing colors for choice. Diameter of head: 46mm for adult use.



### DT113 DUAL HEAD STETHOSCOPE WITH NON-CHILL RING(PEDIATRIC)

Anodized aluminum dual chestpiece matched non-chill rings with a wide range of tubing colors for choice. Diameter of head: 38mm for child use.



### DT114 NEONATAL AND BABY STETHOSCOPE

Zinc alloy chestpiece matched with a wide range of tubing colors for choice.



### DT210 SPRAGUE RAPPAPORT STETHOSCOPE

Chrome-plated brass binaural and full rotatable and various color chestpiece. Accessory kit includes: 2 Diaphragms, 2 Pairs of ear tips, 3 Bells. 22" tube colorful and transparent tube latex free tube for your choice.



**DT211 GOLDEN SPRAGUE RAPPAPORT STETHOSCOPE**  
 Attractive gold-plated head, more deluxe, 22" tube, various color available. Accessory kit includes: 2 Diaphragms, 2 pairs of ear tips, 3 bells.



**DT213 SINGLE TUBE SPRAGUE RAPPAPORT STETHOSCOPE**  
 Chrome-plated brass binaurals with inner spring and rotatable chestpiece.



**DT214 DELUXE SINGLE TUBE SPRAGUE RAPPAPORT STETHOSCOPE**  
 Chrome-plated brass binaural and full rotatable deluxe chestpiece. Special designed double side diaphragm.



**DT215 DELUXE SPRAGUE RAPPAPORT STETHOSCOPE**  
 Chrome-plated brass binaural and full rotatable deluxe chestpiece. Special designed double side diaphragm. 22" tube. Accessory kit includes: 2 Diaphragms, 2 Pairs of ear tips, 3 Bells. Colorful latex free tube for your choice.



**DT216 SPRAGUE RAPPAPORT STETHOSCOPE WITH CLOCK**  
 Chrome-plated brass binaural and full rotatable chestpiece with clock, easy to record heart beat. Colorful latex free tube for your choice.



**DT217 SPRAGUE RAPPAPORT STETHOSCOPE WITH "A" TUBE**  
 Chrome-plated brass binaural with rotatable chestpiece and "A" tube.

# STETHOSCOPE



## DT310 DELUXE DUAL HEAD STETHOSCOPE WITH NON-CHILL RING

Anodized aluminum dual chestpiece matched non-chill rings with a wide range of tubing colors for choice. Inner-spring binaural tubing assembly.



## DT311 DUAL HEAD STETHOSCOPE WITH INNER-SPRING BINAURAL

Anodized aluminum dual head stethoscope with inner-spring binaural, special design suit the ear structure.



## DT312 DELUXE DUAL HEAD STETHOSCOPE WITH NON-CHILL RING

Zinc alloy dual chestpiece matched non-chill rings with a wide range of tubing colors for choice. Inner-spring binaural tubing assembly. More deluxe for professional purpose.



## DT313 SPECIAL SINGLE HEAD STETHOSCOPE

Zinc alloy chestpiece, special design. Comes with inner spring binaural tubing assembly.



## DT314 DELUXE SINGLE HEAD STETHOSCOPE

Zinc-alloy chestpiece with chrome-plated silver or nickel-plated silver finish, laser engraving brand is available. Comes with inner-spring binaural for suiting the ear structure.



## DT315 CROSS-SHAPED MEMBRANE LOCK SINGLE HEAD STETHOSCOPE

Zinc-alloy chestpiece with chrome-plated silver or nickel-plated silver finish. Comes with inner-spring binaural tubing for suiting the ear structure. Special designed cross-shaped membrane lock for getting better heart sound effect through press the back of chestpiece.



**DT316 HUMP SINGLE HEAD STETHOSCOPE**

Zinc alloy hump shape chestpiece. Comes with non-chill ring with inner spring binaural for suiting the ear structure.



**DT317A/B SINGLE HEAD STETHOSCOPE WITH CLOCK**

Chrome-plated brass binaural and anodized aluminum chestpiece with clock, easy to record heart beat.



**DT410A/B CARDIOLOGY STAINLESS STEEL STETHOSCOPE**

Cardiology type, feature a deep cone-shaped chestpiece, comes with non-chill ring. With inner spring binaural tubing assembly. Special designed binaural for suiting the ear structure. With uniquely designed double and thick wall tubing assembly. For adult use.



**DT411 CLASSIC STAINLESS STEEL STETHOSCOPE**

Classic stainless steel chestpiece, with colorful tubing assembly. Special designed inner spring binaural for suiting the ear structure. For adult use.



**DT412 CLASSIC STAINLESS STEEL STETHOSCOPE (PEDIATRIC)**

Classic stainless steel chestpiece, with colorful tubing assembly. Special designed inner spring binaural for suiting the ear structure. For pediatric use.



**DT413 CLASSIC STAINLESS STEEL STETHOSCOPE (INFANT)**

Classic stainless steel chestpiece, with colorful tubing assembly. Special designed inner spring binaural for suiting the ear structure. For infant use.

# STETHOSCOPE



**DT414 STAINLESS STEEL STETHOSCOPE (NEONATE)**  
 Stainless steel chestpiece, with colorful tubing assembly. Special designed inner spring binaural for suiting the ear structure. For neonate use.



**DT415 DOUBLE DIAPHRAGM CARDIOLOGY STAINLESS STEEL STETHOSCOPE**  
 Cardiology type, feature a deep cone-shaped chestpiece, with double side diaphragm and comes with non-chill ring. Special design of binaural for suiting the ear structure. With uniquely designed double and thick wall tubing assembly.



**DT416 STAINLESS STEEL SINGLE HEAD STETHOSCOPE**  
 Deluxe stainless steel chestpiece of diameter 46mm for adult use. Special designed inner spring binaural for suiting the ear structure.



**DT417 CROSS-SHAPED MEMBRANE LOCK STAINLESS STEEL STETHOSCOPE**  
 Super deluxe, special designed cross-shaped membrane lock for getting better heart sound effect through press the back of chestpiece. Comes with inner-spring binaural tubing for suiting the ear structure.



**DT418 THREE PARTS STAINLESS STEEL STETHOSCOPE**  
 Classic stainless steel stethoscope with 3 kinds of chestpiece (adult, pediatric, infant). Comes with inner-spring binaural tubing for suiting the ear structure.



**DT510 ACRYLIC STETHOSCOPE**  
 Acrylic material chestpiece, transparent and deluxe appearance. with inner-spring stainless steel or chrome-plated brass binaural and soft ear tips for suiting the ear structure.






**DT511A/B TEACHING STETHOSCOPE**



**DT512 DISPOSABLE STETHOSCOPE**

Plastic materials, different color for your choice. Especial for disposable use or as a toy.



**DT513A/B OUTDOOR STETHOSCOPE**

Special designed plastic chestpiece provide high quality acoustic effect even in a noisy environment, comes with inner-spring binaural.



**DT514 FETAL STETHOSCOPE**

Special designed chestpiece for pregnancy use. Both getting the heart sound of fetal and pregnant woman by rotate chestpiece.



**Different Packing**

# STETHOSCOPE & NEUROLOGICAL HAMMER



## STETHOSCOPE PARTS

Ear Tip, Name Tag, Tube, Diaphragm, etc.

## NEUROLOGICAL HAMMERS



DT515A1



DT515A2



DT515B



DT515E



DT515D1



DT515D2



DT515F



DT515C



DT515G

## FOUR OF THE MOST POPULAR HAMMERS FOR DIAGNOSIS OF NEUROLOGICAL FUNCTION.

### (DT515A) TAYLOR HAMMER FEATURES:

- The most popular neuro hammer
- Chrome handle or plastic-coating handle with triangular synthetic rubber head in choice of colors
- Latex-free

### (DT515B) BUCKY HAMMER FEATURES:

- Chrome-plated handle with 2 sided synthetic rubber head
- Handle conceals screw-in brush
- Needle contained within the head
- Latex-free

### (DT515D) BABINSKI HAMMER FEATURES :

- Rubber head may be used in both the horizontal or vertical positions
- Latex-free

### (DT515E) WARTENBERG PINWHEEL FEATURES :

- Stainless steel handle with textured grip
- Rotating spur
- Latex-free



HAMMER SET A



HAMMER SET B



HAMMER SET C

# SPHYGMOMANOMETER



## SW1211 PROFESSIONAL MERCURIAL SPHYGMOMANOMETER

Desk type 300mm Hg calibrated scales, cotton cuff, Latex inflation system. Metal case with textured backed enamel finish

|         |                |
|---------|----------------|
| SW1211A | Colour: Gray   |
| SW1211B | Colour: Blue   |
| SW1211C | Colour: Green  |
| SW1211D | Colour: Yellow |
| SW1211E | Colour: Red    |
| SW1211F | Colour: Black  |



**SW1212  
WALL TYPE MERCURIAL  
SPHYGMOMANOMETER**



**SW1213  
FLOOR TYPE MERCURIAL  
SPHYGMOMANOMETER**



## SW1221 ANEROID SPHYGMOMANOMETER (STOP PIN OR NO STOP PIN)

300mm Hg gauge, cotton or nylon cuff with artery label and gauge holder, latex inflation system, complete with leatherette zippered carrying case and operating instruction.



## SW1223 ANEROID SPHYGMOMANOMETER FOR CHILD

200mm Hg gauge, cotton or nylon cuff with artery label and gauge holder, latex inflation system, complete with leatherette zippered carrying case and operating instruction.

# SPHYGMOMANOMETER



SW1224 COLOUR ANEROID SPHYGMOMANOMETER



SW1225 ANEROID SPHYGMOMANOMETER WITH SEPARATE OR ATTACHED STETHOSCOPE



SW1226 ARM TYPE ANEROID SPHYGMOMANOMETER



SW1227 PLASTIC GAUGE ANEROID SPHYGMOMANOMETER



SW1231 PALM TYPE ANEROID SPHYGMOMANOMETER  
Single or double tube available  
Small gauze or Big gauge available



SW1232 PALM TYPE ANEROID SPHYGMOMANOMETER  
Palm type aneroid sphygmomanometer with plastic ring  
Single or double tube available

# SPHYGMOMANOMETER



**SW1237 DESK TYPE ANEROID SPHYGMOMANOMETER**



**SW1237 DESK TYPE ANEROID SPHYGMOMANOMETER**

Fixed or Turning gauge available



**SW1239 FLOOR TYPE ANEROID SPHYGMOMANOMETER**

Round or square gauge available



**SW1238 WALL TYPE ANEROID SPHYGMOMANOMETER**



**SW1242 ARM-TYPE FULLY AUTOMATIC ELECTRONIC BLOOD PRESSURE MONITOR**

Measuring range: 2.7kpa~40kpa (20mmHg~300mmHg)  
Accuracy: 0.5kPa (4mmHg)  
Power source: Type AA batteries x4  
Pressurization: By an air pump (imported pump)  
Deflation: Automatic exhaust  
Automatic power off: To be automatically cut off after 5 minutes of nonuse to save energy.




**SW1241 WRIST-TYPE FULLY AUTOMATIC ELECTRONIC BLOOD PRESSURE MONITOR**

Measuring range: 2.7kpa~40kpa (20mmHg~300mmHg)  
Accuracy: 0.5kPa (4mmHg)  
Power source: Type AAA batteries x2  
Pressurization: By an air pump (imported pump)  
Deflation: Automatic exhaust  
Automatic power off: To be automatically cut off after 5 minutes of nonuse to save energy.


# ECG PRODUCTS



PF, with gel pocket



PF, flat pad



PF, flat pad, ellipse

## SW1402A DISPOSABLE ECG ELECTRODE.

### POLYURETHANE FOAM(PF)

Thickness 1mm.

Adult: Dia 55mm, 50mm


Child: Round Dia 45mm, 40mm

Pyriform: 35x30mm, 48x40mm


Rectangle: 49x47mm

With gel pocket or flat pad

With or without skin abrasive




NW, with gel pocket A




NW, with gel pocket B



NW, pyriform



NW, flat pad



NW, rectangle

## SW1402B DISPOSABLE ECG ELECTRODE. NON-WOVEN(NW)

Adult: Dia 55mm, 50mm

Child: Round Dia 45mm, 40mm, Pyriform 48x40mm, Rectangle 48x47mm

With gel pocket or flat pad

With or without skin abrasive





Disposable Electrode of Cerebro Vascular Function Treatment, Polyurethane Foam, thickness 1mm, 42x24mm. With flat pad.



Magnetic Resonance (MR) quatrode electrodes. Dia 55mm, etc.



Disposable neutral electrode It is a negative plate of the high frequency electroto. It is used for electric welding, electric cutting, etc. of the high frequency electroto. 16X24cm. Single or double pole



Active Carbon Fibre Electrodes, Connector dia 2.0 or 2.2mm 130x90mm, 40x40mm, 50x50mm, 40x100mm, 50x100mm, 60x80mm, 127x76mm, etc. Physical therapy use



### SW1601 ENT DIAGNOSTIC SET, WG-II

Standard set:

- Ophthalmoscope-----1
- Otoscope Carrier-----1
- Ear Speculum(3, 4, 5 and 6mm, diam)-----4
- Nasal Speculum(8mm, diam)-----1
- Dilating Nasal Speculum-----1
- Laryngeal Mirror Stem-----1
- Laryngeal Mirrors-----3
- Tongue Depressor, Plastic-----1
- Battrey Handle-----1
- Spare Bulbs(2large,4small)-----6



### DIAGNOSTIC SET, WG-XP

Standard set:

- Ophthalmoscope-----1
- Otoscope Carrier-----1
- Ear Speculum(2.5, 3, 4, 5mm, diam)-----4
- Dilating Nasal Speculum-----1
- Laryngeal Mirror Stem-----1
- Laryngeal Mirrors-----2
- Tongue Depressor, Plastic-----1
- Battrey Handle-----1
- Spare Bulbs(2.5V/0.3A)-----5



### DIAGNOSTIC SET, WG-III


Standard set:

- Ophthalmoscope-----1
- Otoscope Carrier-----1
- Ear Speculum(2, 3, 4mm, diam)-----3
- Nasal Speculum(8mm, diam)-----1
- Dilating Nasal Speculum-----1
- Laryngeal Mirror Stem-----1
- Laryngeal Mirrors-----3
- Tongue Depressor, Plastic-----1
- Battrey Handle-----1
- Spare Bulbs(2.5V/0.3A)-----4




### SW1401/101 ECG MACHINE

- Input Impedance  $\geq 50M \Omega$
- Input Circuit Current  $\leq 50nA$
- Noise Level  $\leq 15 \mu V_{p-p}$
- Patient Current Leakage  $\leq 10 \mu A$
- CMRR  $\geq 100dB$
- Frequency Response 0.05~150Hz(-3dB)
- A/D Conversion 12bit
- Time Constant >3.2 second
- Skin Voltage  $\pm 300mV$
- Sensitivity 5mm/mV, 10mm/mV  
20mm/mV, Auto
- Paper Speed 25mm/s, 50mm/s
- Power Requirements AC: 220/110V; 50/60HZ  
DC: Rechargeable battery  
supporting 50 exams  
1, 3, 6, 12 channels available



### SW1403 ECG PAPER

50x16x20m  
50x16x30m



### SW1405 THERMAL PAPER

45x50mm  
57x50mm  
80x50mm  
110x50mm , etc.

# THERMOMETER

## CLINICAL THERMOMETER (GLASS)

| ITEM   | DESCRIPTION             |
|--------|-------------------------|
| SW1511 | Rectal Use, short bulb  |
| SW1512 | Oral use, long bulb     |
| SW1513 | Armpit use, small size  |
| SW1514 | Armpit use, middle size |
| SW1515 | Armpit use, large size  |

Back available with: yellow, white, blue  
Scale: °C, F or °C&F



## DIGITAL THERMOMETER

| ITEM    | DESCRIPTION             |
|---------|-------------------------|
| SW1521C | Pen-type, 3-digital, °C |
| SW1521F | Pen-type, 3-digital, F  |
| SW1522C | Pen-type, 4-digital, °C |
| SW1522F | Pen-type, 4-digital, F  |
| SW1523C | Water resistance, °C    |
| SW1523F | Water resistance, F     |
| SW1524C | Flexible tip, °C        |
| SW1524F | Flexible tip, F         |
| SW1528C | Infant use, °C          |
| SW1528F | Infant use, F           |



**SW1532**  
WET DIGITAL THERMOMETER




**SW1531**  
INFRA-RED EAR THERMOMETER



# GYNAECOLOGY & PEDIATRICS



**SW2101A1 DISPOSABLE VAGINAL SPECULUM LATERAL SCREWTYPE**  
**SW2101A2 DISPOSABLE VAGINAL SPECULUM CENTRAL SCREWTYPE**  
 Size: S, M, L



**SW2101B DISPOSABLE VAGINAL SPECULUM PUSH TYPE**  
 Size: S, M, L




**SW2101C1, SW2101C2 FASTEN TYPE**  
 Size: XS, S, M, L

**SW2101D TACHE TYPE**



**SW2102 CERVICAL SCRAPER (WOODEN)**  
**SW2102A WITH ROUND TIP**  
**SW2102B WITH SHARP TIP**




**SW2109 PAP SMEAR KIT**




**SW2103 DISPOSABLE AMNIOTIC MEMBRANE HOOK**



**SW2104 UMBILICAL CORD CLAMP**  
 Color: White, Blue



**SW2105 CERVICAL BRUSH**



**SW2106 INTRAUTERINE CONTRACEPTIVE DEVICE (IUD)**



**SW2108 DISPOSABLE UNDER PAD**  
 Size: 40x60cm, 60x60cm, 60x90cm, 75x75cm, etc



# DRESSING


# PLASTER




**SW3101, SW3102 ZINC OXIDE PLASTER**  
Width: 1.25cm, 2.5cm, 5cm, 7.5cm, 10cm, etc.  
Length: 1m, 5m, 10m, etc.  
Color: white, flesh  
Simple packing, plastic shell available, tinplate shell available



**SW3103 TRANSPARENT SURGICAL TAPE**  
Width: 1.25cm, 2.5cm, 5cm, 7.5cm, etc.  
Length: 5Y, 10Y, etc.




**SW3104 MICROPOROUS SURGICAL TAPE**  
Width: 1.25cm, 2.5cm, 5cm, 7.5cm, etc.  
Length: 5Y, 10Y, etc.  
With or without dispenser



**SW3105 SILK TAPE**  
Width: 1.25cm, 2.5cm, 5cm, 7.5cm, etc.  
Length: 5Y, 10Y, etc.



**SW3106 FIRST-AID PLASTER**  
Size: 72x19mm, etc.  
Color: flesh  
Elastic fabric, PVC, etc.



**SW3108 SELF-ADHESIVE WOUND DRESSING**  
Size: 5x7cm, 9x10cm, 10x15cm, etc.



**SW3112 WOUND DRESSING STRIP**  
Width: 4cm, 6cm, 8cm, etc.  
Length: 1m, 5m, etc.  
With or without hole  
Color: flesh



**SW3115 ELASTOPLASTER**  
Width: 5cm, 7.5cm, 10cm, etc.  
Length: 5m, 10m, etc.  
Color: flesh



## SW3201 LAPAROTOMY SPONGE

40's x40's cotton yarn  
 With or without X-Ray thread  
 With or without blue tape  
 Mesh: 19x15, 24x20, etc.  
 Size: 30x30cm, 45x45cm, etc.  
 Ply: 4, 6, 8, 12, etc.



## SW3202 GAUZE SWAB

40's x40's cotton yarn  
 With or without X-Ray thread  
 Mesh: 19x15, 24x20, etc.  
 Size: 5x5cm, 7.5x7.5cm, 10x10cm, 10x20cm, etc.  
 Ply: 8, 12, 16, etc.  
 Color: White, Green



## SW3203 GAUZE SWAB

(Individual pack)  
 40's x 40's cotton yarn  
 Sterile or non-sterile  
 With or without X-Ray thread  
 Mesh: 19 x 15, 24 x 20, etc.  
 Size: 5x5cm, 7.5x7.5cm, 10x10cm, 10x20cm, etc.  
 Ply: 8, 12, 16, etc.  
 Color: White, Green  
 Pack: Peel pack




## SW3204 GAUZE ROLL

40's x40's cotton yarn  
 With or without X-Ray thread  
 Mesh: 19x15, 24x20, etc.  
 Size: 36" x 50M, 36"x100 Yards, 4Ply etc.



## SW3205 COTTON BALL

Weight: 0.5g/pc, 1g/pc



## SW3206 WOW GAUZE BANDAGE

40's x40's cotton  
 Mesh: 19x15, 24x20, etc.  
 Size: 5cm x4m, 7.5cm x4m, 10cm x4m, etc.

# COTTON PRODUCTS



## SW3207 GAUZE SWAB FOR TRACHEOTOMY

Size: 2"x2", 3"x3", 4"x4"  
 Ply: 4, 6, 12, etc  
 Material: gauze, non-woven



## SW3208 COTTON ROLL

Spec.: 50g, 100g, 200g, 250, 300g, 400g, 454g, 500g, 1000g, etc per roll



## SW3211 DENTAL SPONGES

Diameter: 0.8cm, 1cm, 1.2cm, 1.4cm, etc  
 Length: 3.8cm



## SW3212 SURGICAL TOWEL

Size: 17"x27"  
 Blue, Green, White color  
 with hole 1", 2", 3", 4", etc.



## SW3215 BASIC DRESSING SET

To clean wound and change drug for patients who is suffered from trauma.



## SW3209 ALCOHOL PREP PADS

Contents: 70% isopropyl alcohol  
 Non-woven size: 65x30mm  
 Pouch size: 50x50mm



## SW3216 TRACHEOSTOMY CARE KIT



## SW3217 SUTURE REMOVAL KIT



### SW3301 HIGH ELASTIC BANDAGE

High elastic force tabby bandage  
Common width: 2", 3", 4", 6", 8", etc.  
Stretch rate: 1:2.9~3.1  
Color: natural white and complexion  
Weight: 75g/m<sup>2</sup>, 90g/m<sup>2</sup>



### SW3302 ALL COTTON WRINKLE ELASTIC BANDAGE

Common width: 2", 3", 4", 6", 8"  
Common weight: 75g/m<sup>2</sup>, 90g/m<sup>2</sup>, 100g/m<sup>2</sup>, 140g/m<sup>2</sup>  
Stretch rate: ≥200%



### SW3303 SPANDEX ELASTIC BANDAGE

Common width: 2", 3", 4", 6"  
Stretch rate: ≥200%  
Color: natural white, flesh, bleached or not crepe or plain  
Weight: 75~80g/m<sup>2</sup>



### SW3304 CONFIRMING (PBT) ELASTIC BANDAGE

Common width: 2", 3", 4", 6", etc.  
Common weight: 32±5g/m<sup>2</sup>, thick type 56g/m<sup>2</sup>  
Stretch rate: ≥200%  
Color: natural white, complexion/flesh, bleached



### SW3305 SELVAGE GAUZE BANDAGE

Width: 5cm, 7.5cm, 10cm, 15cm, etc.  
Length: 5m, etc.  
Mesh: 30x20



### SW3306A PLASTER OF PARIS BANDAGE

**SW3306B FIBRE ORTHOPAEDIC BANDAGE**  
Width: 5cm, 7.5cm, 10cm, 15cm, 20cm, etc.  
Length: 2.7m, 3.6m, 4.6m, etc.



### SW3307 SELF-ADHESIVE ELASTIC BANDAGE

Width: 5cm, 7.5cm, 10cm, 15cm, etc.  
Length: 4.5m etc.  
Material: cotton, non-woven



### SW3308 TRIANGULAR BANDAGE

Size: 36"x36"x51", 40"x40"x56", etc.  
Mesh: 40x34, 50x30, 48x48, etc.

# NON-WOVEN PRODUCTS



## SW3411 BOUFFANT CAP

Size: 18", 21", 24", etc.  
Color: Green, Blue, White, etc.



## SW3412 DOCTOR CAP

1ply or 2ply available  
Color: Green, Blue, White, etc.



## SW3413 STRIP CAP

Size: 18", 21", 24", etc.  
Color: Green, Blue, White, Red, etc.




## SW3451 SHOE COVER

Material: PE, CPE, PP, etc.  
Size: 36x15cm, 38x16cm,  
40x17cm, etc.  
Non-skid available



## SW3421 NON-WOVEN FACE MASK

2ply, 3ply or with shield available  
With earloops or ties  
Color: Green, Blue, White, Pink, etc.



## SW3422 DUST-PROOF FACE MASK

Prevents spread of germs and filters out dust and pollen  
Color: White, Blue, Green



## SW3423 PAPER FACE MASK

1ply or 2ply available



## SW3431 SURGICAL GOWN

Material: PP, SMS, Spunlace, etc.  
Size: S, M, L, XL



## SW3415 SLEEVE COVER

Material: PP, PE, etc.



## SW3432 COVERALL

Material: PP, SMS, Spunlace, etc.  
Size: S, M, L, XL



# HEALTH-CARE



# SURGICAL PRODUCTS



## SW4101 SURGICAL BLADE

SW4101A Carbon steel

SW4101B Stainless steel

Size: 10, 11, 12, 12B, 13, 14, 15, 15C, 16, 17, 18, 19, 20, 21, 22, 23, 24, 25, 36#




## SW4102 SCALPEL (PLASTIC HANDLE)

SW4102A Carbon steel

SW4102B Stainless steel

Size: 10, 11, 12, 12B, 13, 14, 15, 15C, 16, 17, 18, 19, 20, 21, 22, 23, 24, 25, 36#



## SW4103 STITCH CUT BLADE

Size: 5cm, 10cm

## SW4104 SKIN GRAFT BLADE

Sterilization With Gamma Radiation



## SW4105 SURGICAL SUTURES

Material: Silk, Nylon, Polypropylene, Catgut, PGA, etc.



## SW4106 TOURNIQUET

Color: orange, blue, green, etc.

Latex tourniquet available



## SW4109 FIRST AID KIT




## SW4110 OXYGEN REGULATOR



## SW4111 PEAK VELOCITY METER



**SW4210 TRACHEOSTOMY MASK**



**SW4211 MEDIUM CONCENTRATION OXYGEN MASK**  
 Size: Adult, adult elongated  
 paediatric, paediatric elongated



**SW4212 NEBULIZER OXYGEN MASK**  
 Size: Adult, adult elongated  
 paediatric, paediatric elongated  
 Volume of nebulizer: 6ml, 15ml, 20ml



**SW4213 NEBULIZER WITH MOUTH PIECE**  
 Volume: 6ml, 15ml, 20ml



**SW4214 HIGH CONCENTRATION OXYGEN MASK**  
 Size: Adult, adult elongated  
 paediatric, paediatric elongated



**SW4215 MULTI-VENT MASKS**  
 Green for 24%, 26%, 28% and 30%  
 White for 35%, 40% and 50%



**SW4216 INHALATOR**  
 Inhaler volume: 180cc / 320cc  
 Size: Infant, Child, Adult  
 Protect from the static electricity can be chosen




**SW4217 CPR MASK**

# TUBE & BAG



**SW4218 VENTURI OXYGEN MASK**  
 Size: Adult, adult elongated  
 paediatric, paediatric elongated  
 with 7 dilutors



**SW4219 NASAL OXYGEN CANNULA**  
 Tube: 7feets with joint  
 Size: XS, S, M, L




**SW4220 COLOR NASAL OXYGEN CANNULA**  
 Tube 50cm



**SW4221 ENDOTRACHEAL TUBE**  
 Model: standard, reinforced,  
 preformed oral, preformed nasal.  
 With cuff: I.D. 3.0-11mm  
 Without cuff: I.D. 2.0-11mm




**SW4229 OXYGEN CANNULA**  
 Tube 45cm



**SW4230 LARYNGEAL MASK DEVICE**



**SW4244 INFANT MEDICINE FEEDER**



**SW4245 INTUBATING STYLET**




**SW4246 TRACHEOSTOMY TUBE**  
 With cuff: I.D. 5.0~10.0mm  
 Without cuff: I.D. 3.0~10.0mm



**SW4222 AIRWAY**

| SIZE       | COLOUR |
|------------|--------|
| 000#, 40mm | pink   |
| 00#, 50mm  | blue   |
| 0#, 60mm   | black  |
| 1#, 70mm   | white  |
| 2#, 80mm   | green  |
| 3#, 90mm   | yellow |
| 4#, 100mm  | red    |
| 6#, 120mm  | orange |



**SW4223 SUCTION CATHETER**  
 Size: 5 Fr~24 Fr  
 PE or Peel packing  
 With or without X-Ray detectable thread



**SW4224 STOMACH TUBE**  
 Size: 6 Fr~24 Fr  
 PE or Peel packing  
 With or without X-Ray detectable thread



**SW4225 FEEDING TUBE**  
 Size: 4 Fr~22 Fr  
 PE or Peel packing  
 With or without X-Ray detectable thread



**SW4226 NELATON TUBE**  
 Size: 6 Fr~24 Fr  
 Female catheter also available




**SW4227 RECTAL TUBE**  
 Size: 24 Fr~36 Fr  
 PE or peel packing




**SW4228 THORACIC CATHETER**

# TUBE & BAG



### SW4231 FOLEY CATHETER

One-way: red, 8 Fr~28 Fr  
 Two-way: 6 Fr~28 Fr  
 Three-way: 16 Fr~28 Fr  
 Soft or hard valve available



### SW4232 ALL SILICONE FOLEY CATHETER

100% silicone material



### SW4235 LATEX TUBING

I.D.: 3mm~12mm  
 O.D.: 5mm~18mm



### SW4236 LATEX T-DRAINAGE TUBE

Sterile, 10~26 Fr



### SW4237 CONDOM CATHETER

SIZE: S, M, L, XL  
 PE or Peel packing



### SW4238 PENROSE DRAINAGE TUBE

Size: 3/8", 1/2", 3/4", 5/8", 1", 2", etc



### SW4241 PAEDIATRIC URINE COLLECTOR

Size: 100ml, 200ml, etc.



### SW4242 DISPOSABLE COLOSTOMY BAG


Size: 20x14cm, 25x15cm



**SW4251 URINE DRAINAGE BAG**  
 Size: 1500ml, 2000ml, etc  
 Specifications: 90/120cm tube, with anti-reflux valve and bottom outlet, Pull-Push valve.



**SW4252 URINE DRAINAGE BAG**  
 Size: 1500ml, 2000ml, etc  
 Specifications: 90/120cm tube, with anti-reflux valve and bottom outlet, Twist valve.



**SW4253 URINE DRAINAGE BAG**  
 Size: 1500ml, 2000ml, etc  
 Specifications: 90/120cm tube, with anti-reflux valve and bottom outlet, T valve.



**SW4250 URINE DRAINAGE BAG**  
 Size: 1500ml, 2000ml, etc  
 Specifications: 90/120cm tube, without bottom outlet



**SW4255 URINE DRAINAGE BAG**  
 Size: 2000ml  
 Specifications: front film transparent, back film white. with needle sample port & plastic tie, anti-reflux valve Twist valve outlet



**SW4257 URINE DRAINAGE BAG**  
 Size: 2000ml  
 Specifications: with anti-reflux valve, needle sample port, reinforced single hanger, finger guard clamp




**SW4259 URINE METER**  
 Size: 2600ml drainage bag, 400ml urine meter  
 Specifications: Easy to read molded plastic meters are calibrated to 400cc. with anti-reflux valve, needle sample port, reinforced double hanger and rope hanger



**SW4261 SUCTION CONNECTING TUBE**  
**SW4262 YANKAUER HANDLE**  
 Plain tip or crown tip



**SW4271 MUCUS EXTRACTOR**



**SW4264 MOUTH PIECE**



**SW4263 WOUND DRAINAGE RESEVOIR**  
 Volume: 200ml, 400ml

# GLOVES



## SW4301 PE GLOVES

High & Low density polyethylene  
Dispenser pack available  
Two finger gloves available



## SW4304 SURGICAL GLOVES

Size: 6#~8.5#  
Powdered and powdered-free available

## SW4305 EXAM GLOVES

Size: S, M, L  
Powdered and powdered-free available



## SW4306 PVC GLOVES

Size: S, M, L  
COLOUR: Yellow, Blue, Green, etc.  
Powdered and powdered-free available




## SW4351 KIDNEY DISHES

Size: 243x125x50mm, 203x102x44mm, 165x80x38mm



## SW4352 STAINLESS STEEL MEDICAL SQUARE TRAY

Size: 300x200mm, 350x250mm, 400x300mm  
450x300mm, 400x350mm, 400x500mm  
400x600mm



## SW4353 STAINLESS STEEL INJECTOR CASE (TROMOL)

Size: 200x80x65mm, 185x80x65mm



## SW4354 STAINLESS STEEL STERILE PLATE

Size: 300x200mm, 12 inch    250x160mm, 9 inch  
200x110mm, 8 inch    150x100mm, 6 inch



## SW4355 STAINLESS STEEL COLLECTING TANK (TROMOL)

Size: 170x140mm, 200x155mm, 230x175mm  
260x190mm, 300x200mm



# INJECTION PRODUCTS



## SW4401 DISPOSABLE SYRINGE

Size: 1ml, 2ml, 3ml, 5ml, 10ml, 20ml, 50ml, etc.  
 Specifications: with or without needle, luer slip or luer lock.  
 3 parts and 2 part available  
 Insulin syringe, Tuberculin syringe available.



## SW4402 TUBERCULIN



## SW4403 INSULIN 100 units



## 100200A 100130A 100150A 200150B Accessory

SW4406/100200A CT MEDRAD MCT PLUS. VISTRON CT. ENVISION  
 SW4406/100130A DSA MEDRAD MARK IV  
 SW4406/100150A DSA MEDRAD MARK V  
 SW4406/200150B DSA LF ANGIOMAT 6000  
 SW4406/10060A MRI MEDRAD SPECTRIS & SOLARIS  
 SW4406/10060C MRI MEDRAD SPCCTRIS & SOLARIS



## SW4404 SAFE TYPE SYRINGE



## SW4407 IRRIGATION SYRINGE



## SW4416 I.V. CANNULA

Size: 14G, 16G, 18G, 20G, 22G, 24G  
 A: with injection valve  
 B: Butterfly with wings  
 C: Pen like



## SW4417 3WAY STOPCOCK

With or without extension tube




## SW4419 HEPARIN CAP

# INJECTION PRODUCTS



### SW4408 DISPOSABLE INFUSION SET

Specifications: vented spike with cap, plastic chamber, PVC tube 150cm, with or without needle.




### SW4409 DISPOSABLE BLOOD TRANSFUSION SET

Specifications: vented spike with cap, PVC tube 150cm, with or without hypodermic needle.



### SW4410 CHILDREN HANGING TRANSFUSION DEVICE

Size: 100ml, 150ml, etc.



### SW4411 DISPOSABLE SCALP VEIN SET

Specifications: butterfly needle 18G~27G, luer slip or luer lock, poly-bag or blister packing.



### SW4412 BLOOD BAG

Size: 250ml, 450ml, etc.  
Single, Double, Triple, Quadruple bag available.



### SW4413 DISPOSABLE DENTAL NEEDLE

| Specification               | 27G(0.4mm)             | 30G(0.3mm)         |
|-----------------------------|------------------------|--------------------|
| Color                       | Brown                  | Yellow             |
| Length(mm)                  | 13, 21, 25, 32, 38, 42 | 13, 21, 25, 32, 38 |
| Hub: Metric System ISO 7885 |                        |                    |
|                             | British System         |                    |



### SW4415 DISPOSABLE HYPODERMIC NEEDLE

Needle size: 18G-27G  
Blister packing or poly-bag.



### SW4418 DISPOSABLE VACUTAINER NEEDLE

Needle size: 20G, 21G, 22G, 23G  
Length: 1.5"



### SW4501 HOT WATER BAG

With or without cover  
Size: 1000ml, 1500ml, 2000ml, etc.



### SW4502 ICE BAG

Size: 4"~11"  
Cloth or rubber type



### SW4503 AIR CUSHION

Size:  $\phi$ 40cm,  $\phi$ 42.5cm,  $\phi$ 45cm



### SW4504 PATIENT ID BRACELETS

Adult, Infant  
Color: Pink, Blue, White, etc  
Insert type or write on type



### SW4506 EAR SYRINGE BALL

Size: 30ml, 60ml, 90ml, etc



### SW4509 INSTANT ICE PACK



### SW4510 HOT/COLD PACK



### SW4505 LATEX CONDOM

| SW4505A        | SW4505B        | SW4505C        |
|----------------|----------------|----------------|
| Smooth surface | Dotted surface | Ribbed surface |
| Lubricated     | Lubricated     | lubricated     |


Size: 49mm, 52mm, 55mm  
Be available with various fruit smell  
Available with various packing  
Specification with EN600, ISO4074, ASTM, WHO, SABS, etc.

# SCALES



## SW4611 BATHROOM SCALE

Capacity: 120kg / 260 lb  
Graduation: 1kg/2 lb



## SW4612 ELECTRONIC BATHROOM SCALE

Capacity: 2~150kg/4~330 lb  
Division: d=0.1kg/0.2 lb  
Indication: 4's Place in LCD of 25mm(1 inch)  
Power: 1\*CR2032 Li-ion battery  
Low Power warning: "Lo"  
Over load warning and automatic "0" place tracing  
Vibrate to turn on and Automatic turn off in 6~8 seconds




## SW4621 BABY BALANCE

Max weight: 10kg  
Min value per division: 50g  
Load platform area(L x W): 545x270mm




## SW4622A/B BABY BALANCE

Max weight: 20kg  
Min value per division: 0~10kg: 50g  
10~20kg: 100g  
Load platform area (L x W): 545x270mm



## SW4625A/B ELECTRONICBABY BALANCE

Max weight: 20kg  
Min value per division: 10g  
Load platform area (L x W): 545x270mm




## SYE-2004E2 KITCHEN SCALE

Max Capacity: 2000g/70.55 oz (3000g is also available)  
Div.: 1g/0.05 oz



### TABLE BALANCE

| ITEMS   | CAPACITY   | SENSITIVITY |
|---------|------------|-------------|
| JPT-1   | Max: 100g  | D: 0.1g     |
| JPT-2   | Max: 200g  | D: 0.2g     |
| JPT-5   | Max: 500g  | D: 0.5g     |
| JPT-10  | Max: 1000g | D: 1g       |
| JPT-20A | Max: 2kg   | D: 2g       |
| JPT-50A | Max: 5kg   | D: 5g       |



### SW4630 ELECTRONICAL TABLE BALANCE

| MAXIUM VOLUME                                      | SENSITIVITY(LED Display) |
|--|--------------------------|
| 100g, 200g, 300g, 400g,<br>500g, 600g,1000g, 2000g | 0.01g                    |
| 600g, 1500g  | 0.1g                     |



### SW4631 WEIGHT HEIGHT BALANCE

Max weight: 120kg  
 Height range: 70~190cm  
 Min weight division: 0.5kg  
 Min value of height division: 5mm  
 Load platform area (L x W): 375x275mm



### SW4632 ELECTEONIC BODY SCALE

Max weight: 200kg  
 Min weight division: 0.1kg  
 Height rang: 70~190cm  
 Min height per division: 5mm  
 Load platform area(LxW): 375x275mm



### SW4633 BODY WEIGHT BALANCE

Max weight: 160~200kg  
 Height range: 75~200cm  
 Scale division(d): 100g  
 Min value of height division: 5mm  
 Load platform area (L x W): 375x275mm



### SW4642 ELECTRONIC WHEELCHAIR BALANCE

Max weight: 200kg  
 Min weight division: 0.1kg

# HOSPITAL BED



## SW4701 HOSPITAL FLAT BED

Size: 2000(length) x 900(width) x 500(height) mm  
Paint Coated type, Stainless steel type



## SW4706 HOSPITAL FOLDING BED WITH THREE PARTS

Size: 2000(length) x 900(width) x 500(height) mm  
Paint Coated type, Stainless steel type



## SW4711 MANUAL 3-LEVER BED WITH WHEELS, FLAT PIPE ALUMINUM PLATE BED HEAD, AND ABS BED SURFACE

- Size: 2000(length) x 900(width) x 500(height) mm
- Flat pipe, aluminum bed head and end, flexible assembly.
  - Horizontal aluminum alloy guardrail, safe and reliable. Can be positioned at vertically or horizontally.
  - Whole-braking double-face pedals for 4 wheels.



## SW4712 DOUBLE-LEVER BED WITH WHEELS, ABS BED HEAD, AND STEEL PLATE SURFACE

- ABS bed head and end, flexible assembly.
- Steel guardrail, safe and reliable. Can be positioned at vertically or horizontally.
- Whole-braking double-face pedals for 4 wheels, strong stability and no noise.
- Double-face pressure bearing in the lever, and joined with universal joints, Easily up and down, stable and reliable.




## SW4718 BEDSORE BLOCK AIR MATTRESS BEDSORE REMEDY AIR MATTRESS



Chromeplate    Stainless steel    Aluminum alloy

## SW4730 TRANSFUSION STAND

Product size (LxWxH): 218x57x42cm(Topmost)  
122x57x42cm(Lowest)



## SW4726 CAMPING BED

Product size (LxWxH): 190x57x40cm  
The bearing pressure: ≤159kg



**SW4723 ALUMINUM ALLOY STRETCHER**  
 Max. Unfolded Size(L x W x H): 202x42x7cm  
 min. Unfolded Size(L x W x H): 168x42x7cm  
 Self-weight: 9kg  
 Load Bearing: 120kg



**SW4719 FOLDING STRETCHER**  
 Product size (LxWxH): 221x53x16cm  
 The bearing pressure:  $\leq 159$ kg



**SW4725 SOFT STRETCHER**  
 Size(L x W): (115~180)cm x62cm  
 Net weight: 2kg  
 The bearing pressure:  $\leq 159$ kg



**SW4720 AUTOMATIC LOADING STRETCHER**  
 Size(L x W x H): 190x54x90cm  
 Bearing pressure: 159kg



**SW4721 ALUMINIUM-ALLOY STRETCHER FOR AMBULANCE**  
 Size(L x W x H): 190x54x33cm  
 Bearing pressure: 159Kg



**SW4722 MULTI-PURPOSE ALUMINUM ALLOY EMERGENCY BED**  
 Size(L x W x H): 195x55x90cm  
 Bearing pressure: 159kg

A black and white photograph of a wheelchair, viewed from a side-rear perspective. The wheelchair is positioned in the center-right of the frame. The background is dark and atmospheric, with several bright, out-of-focus light sources creating a bokeh effect. A white grid is overlaid on the image, consisting of a vertical line on the left, a horizontal line across the middle, and a vertical line on the right. A white rectangular box is centered horizontally, containing the word "HEALING" in a serif font.

# HEALING





SW/FS809



SW/FS874A



SW/FS874B



SW/FS874Y

## WHEELCHAIR

| ITEM NUMBER | MEASUREMENTS (cm) |    |        |    |    |    |    |     |    |    | N.W. (KG) |
|-------------|-------------------|----|--------|----|----|----|----|-----|----|----|-----------|
|             |                   |    |        |    |    |    |    |     |    |    |           |
| SW/FS809    | 60(65)            | 24 | 41(46) | 59 | 20 | 50 | 87 | 106 | 40 | 45 | 19        |
| SW/FS874A   | 60(65)            | 28 | 41(46) | 59 | 20 | 47 | 91 | 106 | 40 | 45 | 19        |
| SW/FS874B   | 60(65)            | 28 | 41(46) | 59 | 20 | 47 | 91 | 106 | 40 | 45 | 19        |
| SW/FS874Y   | 60(65)            | 30 | 41(46) | 59 | 20 | 51 | 91 | 106 | 40 | 45 | 19        |

A: DROP BACK HANDLE    B: MAG WHEEL  
Y: HARD CUSHION SEAT    (46)(51): WIDE SEAT

# WHEELCHAIR



SWIFS874L



SWIFS903L



SWIFS874LAH-35



SWIFS907LABH



SWIFS721L-36

## WHEELCHAIR

| ITEM NUMBER    | MEASUREMENTS (cm) |    |        |    |    |    |    |     |    |    | N.W. (KG) |
|----------------|-------------------|----|--------|----|----|----|----|-----|----|----|-----------|
|                |                   |    |        |    |    |    |    |     |    |    |           |
| SWIFS874L-41   | 63                | 28 | 41     | 57 | 15 | 45 | 82 | 100 | 40 | 45 | 13        |
| SWIFS874L-46   | 66                | 28 | 46     | 61 | 20 | 47 | 88 | 106 | 40 | 45 | 13        |
| SWIFS903L      | 58(63)            | 28 | 41(46) | 61 | 15 | 50 | 90 | 120 | 40 | 45 | 14        |
| SWIFS874LAH-35 | 58                | 28 | 35     | 52 | 15 | 45 | 85 | 93  | 36 | 45 | 12        |
| SWIFS907LABH   | 56(61)            | 28 | 41(46) | 40 | 15 | 45 | 90 | 100 | 40 | 45 | 13        |
| SWIFS721L-36   | 58                | 35 | 36     | 59 | 10 | 48 | 70 | 85  | 38 | 27 | 13        |

A: DROP BACK HANDLE    B: MAG WHEEL    L: ALUMINUM  
H: HAND BRAKE    -35: SEAT WIDTH

# WHEELCHAIR



SW/FS904B



SW/FS976AB-43



SW/FS110A ELECTRIC TYPE



SW/FS693



SW/FS895L

## WHEELCHAIR, COMMODE CHAIR

| ITEM NUMBER   | MEASUREMENTS (cm) |    |        |      |      |       |    |     |    |    | N.W. (KG) |
|---------------|-------------------|----|--------|------|------|-------|----|-----|----|----|-----------|
|               |                   |    |        |      |      |       |    |     |    |    |           |
| SW/FS904B     | 58(63)            | 25 | 41(46) | 30   | 20   | 50    | 88 | 98  | 40 | 45 | 18        |
| SW/FS976AB-43 | 62                | 26 | 43     | 20   | 20   | 50    | 94 | 98  | 40 | 45 | 13        |
| SW/FS110A     | 58(62)            | 35 | 41(46) | 30   | 20   | 50    | 91 | 106 | 40 | 45 | 30        |
| SW/FS693      | 55                |    | 44     | 12.5 | 12.5 | 48    | 84 | 70  | 44 | 45 | 18        |
| SW/FS895L     | 60                |    |        |      |      | 41-51 | 70 | 48  | 44 | 21 | 5         |

A: DROP BACK HANDLE    B: MAG WHEEL    L: ALUMINUM  
 C: ELEVATING LEGREST    G: RECLINING HIGH BACK

# WALKERS



SW/FS912L



SW/FS961L



(L) (M) (S)

SW/FS917L



SW/FS913L

## WALKERS

| ITEM NO.      | LENGTH(CM) | WIDTH(CM) | HEIGHT(CM) | N.W.(KG) |
|---------------|------------|-----------|------------|----------|
| SW/FS912L     | 57         | 60        | 74.5-92    | 2.6      |
| SW/FS913L     | 56         | 57        | 74.5-92    | 2.4      |
| SW/FS917L (L) | 57         | 42        | 87-97      | 1.4      |
| SW/FS917L (M) | 57         | 42        | 77-87      | 1.2      |
| SW/FS917L (S) | 57         | 41        | 67-77      | 1.0      |
| SW/FS961L     | 56         | 60        | 74.5-92    | 2.8      |

L: ALUMINUM



SW/F S922



SW/F S921



SW/F S934



SW/F S932



SW/F S926



SW/F S910L



SW/F S930L



SW/F S929L



SW/F S928L



SW/F S923L

## STICKS

| ITEM NUMBER   | HEIGHT(CM) | N.W.(KG) |
|---------------|------------|----------|
| SW/F S922     | 72.5-95    | 0.9      |
| SW/F S921     | 72.5-95    | 0.9      |
| SW/F S934     | 65-87.5    | 0.5      |
| SW/F S926     | 80-102.5   | 0.9      |
| SW/F S932     | 74-96      | 0.8      |
| SW/F S910L    | 89         | 0.3      |
| SW/F S930L    | 68.5-91    | 0.3      |
| SW/F S929L    | 68.5-91    | 0.3      |
| SW/F S928L    | 77.5-100   | 0.3      |
| SW/F S923L    | 87-117.5   | 0.5      |
| SW/F S927L    | 82-92      | 0.3      |
| SW/F S925L(S) | 103-123    | 0.8      |
| SW/F S925L(M) | 116-136    | 0.85     |
| SW/F S925L(L) | 128-148    | 0.9      |
| SW/F S936L    | 77.5-100   | 0.25     |

L: ALUMINUM



SW/F S936L



SW/F S927L




SW/F S925L(M)


The image is a composite of two photographs. The top photograph shows a modern dental chair with a reclined seat and adjustable armrests. The bottom photograph shows a surgical microscope mounted on a stand with drawers. The microscope has a control panel with various buttons and a digital display. The entire image has a dark, monochromatic color scheme with a white border around the text.

# EQUIPMENT

## DS-1 SERIES ELECTRIC OPERATING TABLE FOR CEREBRALIS SURGERY



DS-1




DS-1A


| ITEM  | DS-1                       | DS-1A                      |
|---|----------------------------|----------------------------|
| Length and width for tabletop                         | 1960mmX480mm               | 1960mmX480mm               |
| Min. & Max. height for tabletop                       | ≤500mm, ≥800mm             | ≤500mm, ≥800mm             |
| Max. degree for tabletop forward and backward tilting | forward≥22° backward≥15°   | forward≥22° backward≥15°   |
| Max. degree for tabletop left and right tilting       | left≥15° right≥15°         | left≥15° right≥15°         |
| Max. degree for head rest adjustment                  | up≥45° down≥90°            | up≥45° down≥90°            |
| Max. degree for leg rest adjustment                   | down≥55° up≥5° outward≥90° | down≥55° up≥5° outward≥90° |
| Max. degree for back rest adjustment                  | up≥75° down≥20°            | up≥75° down≥20°            |
| Power   | ~220V, 50HZ, 1.0KW         | ~220V, 50HZ, 1.0KW         |

Kidney bridge can be raised upon 120mm leg section can be fixed as prearrange poses and angles

## DS-3 SERIES ELECTRIC OPERATING TABLE FOR CEREBRALIS SURGERY



DS-3




DS-3A

| ITEM  | DS-3                     | DS-3A                    |
|---|--------------------------|--------------------------|
| Length and width for tabletop                         | 1960mmX480mm             | 1960mmX480mm             |
| Min. & Max. height for tabletop                       | ≤750mm, ≥950mm           | ≤750mm, ≥950mm           |
| Max. degree for tabletop forward and backward tilting | forward≥22° backward≥15° | forward≥22° backward≥15° |
| Max. degree for tabletop left and right tilting       | left≥15° right≥15°       | left≥15° right≥15°       |
| Max. degree for head rest adjustment                  | up≥45° down≥90°          | up≥45° down≥90°          |
| Max. degree for leg rest adjustment                   | down≥90° outward≥90°     | down≥90° outward≥90°     |
| Max. degree for back rest adjustment                  | up≥75° down≥20°          | up≥75° down≥20°          |
| Power   | ~220V, 50HZ, 1.0KW       | ~220V, 50HZ, 1.0KW       |

Kidney bridge can be raised upon 120mm leg section can be fixed as prearrange poses and angles

# OPERATING TABLE

## DS-5 ELECTRIC OPERATING TABLE FOR TRACTING IN ORTHOPEDICS DEPT



DS-5



1006

| ITEM  | DS-5  |
|---|---|
| Length and width for tabletop                         | 1960mmX480mm  |
| Min. & Max. height for tabletop                       | ≤750mm, ≥950mm  |
| Max. degree for tabletop forward and backward tilting | forward≥22° backward≥15°  |
| Max. degree for tabletop trendelenburg                | left≥15° right≥15°  |
| Max. degree for head rest adjustment                  | up≥45° down≥90°   |
| Max. degree for leg rest adjustment                   | down≥55° up≥5° outward≥90°  |
| Max. degree for back rest adjustment                  | up≥75° down≥20°   |
| Power   | ~220V, 50HZ, 1.0kW  |
| ITEM  | 1006 TRACTION DEVICE  |
| Min. height for floor                                 | ≤900mm  |
| Up and down top range                                 | ≥320mm  |
| The stroke of traction arm                            | ≥190mm  |
| Traction Lead the angle                               | up and down turn the turn adjustment 120°<br>The level turn the adjustment 360° |
| Main unit horizontal adjustment                       | ≥120°   |
| The forward and backward stretch adjustment           | ≥1200mm   |
| Height adjustment                                     | ≥300mm  |

1006 Traction device can be chosen

## 3005/DS-2 OPHTHALMOLOGY OPERATING TABLE



3005



DS-2

| ITEM                              | 3005    | DS-2    |
|-----------------------------------|---------|---------|
| Min. height for tabletop to floor | ≤500mm  | ≤500mm  |
| Tabletop up and down range        | ≥300mm  | ≥300mm  |
| Head rest up and down range       | ≥180mm  | ≥180mm  |
| Arm support up and down range     | ≥240mm  | ≥240mm  |
| Tabletop length                   | ≥1900mm | ≥2000mm |
| Tabletop width                    | ≥520mm  | ≥520mm  |



## DS-7 INTEGRATED ELECTRIC OPERATING TABLE



DS-7




DS-4

| ITEM  | DS-7                     | ITEM                             | DS-4               |
|---|--------------------------|----------------------------------|--------------------|
| Length and width for tabletop                         | 1980mmX500mm             | Min.height for tabletop to floor | 600mm              |
| Min.& Max. height for tabletop                        | ≤1100mm, ≥690mm          | Tabletop up and down range       | ≥280mm             |
| Max. degree for tabletop forward and backward tilting | forward≥30° backward≥30° | Tabletop forward tilt            | forward≥10°        |
| Max. degree for tabletop left and right tilting       | left≥18° right≥18°       | Tabletop backward tilt           | backward≥25°       |
| Max. degree for back rest adjustment                  | up≥30° down≥90°          | Backsection up adjusted          | up≥60°             |
| Max. degree for leg rest adjustment                   | down≥80° outward≥90°     | Backsection down adjusted        | down≥10°           |
| Lateral longitudinal shifting of the tabletop         | ≥340mm                   | Tabletop dimensions              | 1240mmX570mm       |
| Max. Up moving of waist section                       | ≥110mm                   | Auxiliary tabletop dimensions    | 600mmX560mm        |
| Power   | ~220V, 50HZ, 1.0KW       | Power                            | ~220V, 50HZ, 1.0KW |

## 3008/308(A)/3008(S) HERD OPERATING UNIVERSAL TABLE



3008(A)



3008



### 3008(S)

Tabletop extrudes when back section turns back with linkage≤135°

| ITEM  | 3008/3008(A)/3008(S)  |
|---|---|
| Min.& Max. height for tabletop                        | 1040mmX800mm/3008(S)1020mmX780mm                            |
| Max. degree for tabletop forward and backward tilting | forward≥25° backward≥25°<br>3008(S)forward≥22° backward≥22° |
| Max. degree for tabletop left and right tilting       | left≥18° right≥18°  |
| Max. degree for back rest adjustment                  | up≥75° down≥10° /Together down≥15° (S)                      |
| Length and width for tabletop                         | 2100mmX480mm/3008(S)2000mmX500mm                            |
| Max. waist rest up to                                 | 120mm/3008(S) Together up≤135°                              |
| Leg rest down adjusted                                | down≥90°<br>3008(S)down≥85° outward≥90° (S)                 |

The cover, tabletop and bass of table of 3008(A) are made of stainless steel.

# OPERATING LAMP



L2000-3E



L734E



L735E

|                                      | L2000-3E                         | L734E                            | L735E                            |
|--------------------------------------|----------------------------------|----------------------------------|----------------------------------|
| Illuminance                          | 60000±500LUX                     | 40000±500LUX                     | 50000±500LUX                     |
| Color Temperature                    | 4500±500K                        | 4000±500K                        | 4000±500K                        |
| Dia. of Light Spot                   | 100~250                          | 100~200                          | 100~200                          |
| Power Voltage                        | 220/24VVAC/DC                    | 220/24VVAC/DC                    | 220/24VVAC/DC                    |
| Rated Voltage of Bulb                | 24V                              | 24V                              | 24V                              |
| Rated Power of Bulb                  | 50W                              | 25W                              | 25W                              |
| Duration of Extended Power Supplying | 2≥HR                             | 3≥HR                             | 2.5≥HR                           |
| Charge Type of Storage Battery       | Quick charge and floating charge | Quick charge and floating charge | Quick charge and floating charge |



L734-II



L734



L751-II




L751

## L734-II COLD LIGHT SHADOWLESS OPERATING LAMP WITH 4 REFLECTORS (Tungsten Halogen Lamp)

|                       |                     |
|-----------------------|---------------------|
| Illuminance           | L734-II ≥50000Lx    |
|                       | L734 ≥40000Lx       |
|                       | 734 ≥24000Lx        |
| Color Temperature     | L734-II 4500±500K   |
|                       | L734, 734 4000±500K |
| Mains Voltage         | ~220V±22V 50Hz±1Hz  |
| Rated Voltage of Bulb | 24V                 |
| Rated Power of Bulb   | 25W                 |

## L751-II COLD LIGHT OPERATING LAMP WITH ONE REFLECTORS (Tungsten Halogen Lamp)

|                       |                    |
|-----------------------|--------------------|
| Illuminance           | L751-II ≥20000Lx   |
|                       | L751 ≥12000Lx      |
| Color Temperature     | L751-II 3500±500K  |
|                       | 751 3000±500K      |
| Mains Voltage         | ~220V±22V 50Hz±1Hz |
| Rated Voltage of Bulb | 24V                |
| Rated Power of Bulb   | L751-II 50W        |
|                       | L751 25W           |



L7/5




L7(W) WALL-MOUNTED OPERATION LAMP

| TECHNICAL PARAMETER                | SR-L7        | SR-L5        |
|------------------------------------|--------------|--------------|
| Light strength (Max.)              | 150,000LUX   | 90,000LUX    |
| Regulating scope of light strength | 100%~50%     | 100%~50%     |
| Bulb power                         | 150W         | 100W         |
| Facular focus depth                | 70cm         | 70cm         |
| Facular diameter                   | 16 to 25cm   | 12to 18cm    |
| Lamp panel diameter                | 700mm        | 520mm        |
| Color temperature                  | 4300K(+/-5%) | 4300K(+/-5%) |
| C.R.I( color revert index )        | 93           | 93           |
| Sucking percent of infrared ray    | 99%          | 99%          |



L2000 6+3-II



L2000 6-II



### L2000 6+3-II COLD LIGHT SHADOWLESS OPERATING LAMP WITH SUPPLEMENTARY LAMP(TUNGSTEN HALOGEN LAMP)

|                                    |   |
|------------------------------------|---|
| IIIUuminance of Main Lamp          | 80000~120000Lx(Stepless Light Adjustment) |
| IIIUuminance of Supplementary Lamp | 60000Lx(Stepless Light Adjustment)        |
| Color Temperature                  | 4500 ±500K                                |
| Mains Voltage                      | ~220V±22V 50Hz ±1Hz                       |
| Rated Voltage of Bulb              | 24V                                       |
| Rated Power of Bulb                | 50W                                       |
| The lowest Height for Installation | 3000mm                                    |

### L2000 6-II COLD LIGHT SHADOWLESS OPERATING LAMP (TUNGSTEN HALOGEN LAMP)

|                                    |   |
|------------------------------------|---|
| IIIUuminance                       | 80000~120000Lx(Stepless Light Adjustment) |
| Color Temperature                  | 4500 ±500K                                |
| Mains Voltage                      | ~220V±22V 50Hz ±1Hz                       |
| Rated Voltage of Bulb              | 24V                                       |
| Rated Power of Bulb                | 50W                                       |
| The lowest Height for Installation | 3000mm                                    |

# OPERATING LAMP



## ZMD-II COLD LIGHT SHADOWLESS OPERATING LAMP WITH SUPPLEMENTARY LAMP(TUNGSTEN HALOGEN LAMP)

|                                    |   |
|------------------------------------|---|
| IIIUuminance of Main Lamp          | ZMD-II 80000~140000Lx(4 Steps Adjustable)<br>ZMD 80000~120000Lx(4 Steps Adjustable) |
| IIIUuminance of Supplementary Lamp | ZMD-II $\geq 60000$ Lx<br>ZMD $\geq 50000$ Lx                                       |
| Color Temperature                  | ZMD-II $4500 \pm 500$ K<br>ZMD $4000 \pm 500$ K                                     |
| Mains Voltage                      | $\sim 220V \pm 22V$ 50Hz $\pm 1$ Hz   |
| Rated Voltage of Bulb              | 24V   |
| Rated Power of Bulb                | 25W   |



## L7412-II COLD LIGHT SHADOWLESS OPERATING LAMP WITH 12 REFLECTORS(TUNGSTEN HALOGEN LAMP)

|                       |   |
|-----------------------|---|
| IIIUuminance          | L7412-II 80000~120000Lx(4 Steps Adjustable)<br>L7412 80000~120000Lx(4 Steps Adjustable)<br>7412 48000~72000Lx(4 Steps Adjustable) |
| Color Temperature     | L7412-II $4500 \pm 500$ K<br>L7412, 7412 $4000 \pm 500$ K   |
| Mains Voltage         | $\sim 220V \pm 22V$ 50Hz $\pm 1$ Hz   |
| Rated Voltage of Bulb | 24V   |
| Rated Power of Bulb   | 25W   |



## L739-II COLD LIGHT SHADOWLESS OPERATING LAMP WITH 9 REFLECTORS(BUILT-IN TYPE)

|                       |  |
|-----------------------|--|
| IIIUuminance          | L739-II $\geq 90000$ Lx<br>739 $\geq 54000$ Lx |
| Color Temperature     | $4000 \pm 500$ K                               |
| Mains Voltage         | $\sim 220V \pm 22V$ 50Hz $\pm 1$ Hz            |
| Rated Voltage of Bulb | 24V  |
| Rated Power of Bulb   | 25W  |



### SW6201 50MA BED SIDE X-RAY UNIT

|                                      |  |                         |
|--------------------------------------|--|-------------------------|
| Max rated capability                 | Fluoroscopy  | 75kV 3mA                |
|                                      | Radiography  | 90kV 50mA               |
| Power supply                         | Voltage  | 180-240V (single phase) |
|                                      | Frequency  | 50Hz                    |
|                                      | Power  | not lower than 3.5KVA   |
|                                      |  | ≤1 ohm (Resistance)     |
| Fluorescent screen size              | 280mm x 356mm  |                         |
| Range of the timer                   | 0.08~6.3s(divided into steps)                        |                         |
| Distance from focus to screen        | 670mm  |                         |
| Specification of X-ray tube          | Model: XD3-3/100<br>(Stationary anode, single focus) |                         |
| Tube head carriage movement range    | 950mm(Vertical), 160mm(Horizontal)                   |                         |
| X-ray tube focus transverse rotation | ±60°   |                         |


### SW6202 100mA MOBILE X-RAY UNIT

|                                   |                               |                     |     |    |     |
|-----------------------------------|-------------------------------|---------------------|-----|----|-----|
| Max rated capability              | Current mA                    | 15                  | 30  | 60 | 100 |
|                                   | Voltage KV                    | 90                  | 90  | 90 | 80  |
|                                   | Time sec.                     | 6.3                 | 6.3 | 4  | 3.2 |
| Power supply                      | Voltage                       | 180~240V            |     |    |     |
|                                   | Frequency                     | 50Hz                |     |    |     |
|                                   | Power                         | not lower than 4KVA |     |    |     |
|                                   |                               | ≤1 ohm (Resistance) |     |    |     |
| Range of the timer                | 0.08~6.3s(divided into steps) |                     |     |    |     |
| Distance between focus and ground | 680mm~1760mm                  |                     |     |    |     |
| Specification of X-ray tube       | Model: XD4-2.9/100            |                     |     |    |     |
| The X-ray tube                    | Longitudinal movement         | 180mm               |     |    |     |
|                                   | Transverse rotation           | ±60°                |     |    |     |



### X-RAY CASSETTES

Size: 5"x7", 8"x10", 10"x12", 11"x14", 12"x15", 14"x14", 14"x17"



### FILM HANGER


Size: 5"x7", 8"x10", 10"x12", 11"x14", 12"x15", 14"x14", 14"x17"



### INTENSIFIER SCREEN

Blue or green sensitive  
Size: 5"x7", 8"x10", 10"x12", 11"x14", 12"x15", 14"x14", 14"x17"

# X-RAY UNIT



## YZ-300 MEDICAL X-RAY UNIT

|                                     |   |  |
|-------------------------------------|---|--|
| Demand of power supply              | Capacity  | Two-phase 30KVA  |
|                                     | Voltage frequency   | 380V / 220V $\pm$ 10% 50HZ $\pm$ 0.5HZ                                       |
|                                     | Internal resistance   | 0.75 $\Omega$ at 380V; 0.25 $\Omega$ at 220V                                 |
| Power supply                        | Adjustable range  | 380V $\pm$ 10% at 380V; 220V $\pm$ 10% at 220V being adjustable continuously |
|                                     | Photograph  |  |
| Photograph                          | Voltage   | 50~100KV, *50~125KV being adjustable continuously                            |
|                                     | Current   | small focus: 50mA, 100mA big focus: 50mA, 100mA, 200mA, 300mA                |
|                                     | Time  | 0.04~6.3s, 23 grades together grading according to R10 coefficient           |
| High-voltage generator              | Capacity  | 20KV A(instantaneous)  |
|                                     | Max. DC output voltage  | 100KV 125KV  |
|                                     | Max. DC output current  | 200mA 300mA  |
| X-Ray tube                          | Model   | XD51-20, 40/100 XD51-20, 40/125  |
|                                     | Focus   | small focus: 1x1mm; big focus: 2x2mm   |
| Photographic bed(YZ-SCI)            | Length x Width x Height of bed face                           | 2000x710x700mm   |
|                                     | Moving range of bed face                                      | 680mm in length; 200mm in breadth  |
|                                     | Ray-filter of bed   | Travel in length $\geq$ 500mm  |
|                                     |   | Grid density N28   |
|                                     |   | Convergence distance fo100cm<br>Grid ratio r8                                |
| Pillar of X-Ray tube unit           | Moving in length along photographic bed                       | 1800mm or 1600mm   |
|                                     | Distance to ground as moving up and down (along upright post) | 650~1800mm   |
|                                     | Rotary around the center or cross arm                         | $\pm$ 180°   |
|                                     | Rotary around the shaft line of X-Ray annular tubes           | -10° ~60° ~+ 120°  |
| Max size of cassette for photograph | 356x432mm(14" x 17")  |  |

Note: The parameters with \*mark was for model Yz300I 300mA 125KV.



**SW/CS300A  
CHAIR MOUNTED DENTAL UNIT**

Power source: 220V 50Hz  
 Air source: 0.55MPa  
 Water source: 0.2~0.4MPa  
 Gross weight: 165Kg+230Kg  
 Dimensions: 114×102×134(cm), 130×79×123(cm)  
 Can be equipped with:  
 High-speed turbine handpiece  
 Low speed air motor handpiece  
 3-way syringe  
 Operating light  
 Air vacuum  
 Saliva ejector  
 X-ray viewer  
 Dental chair  
 Dental stool  
 Built-in ultrasonic scaler



**SW/CS35C OIL PUMP DENTAL CHAIR**

Cushion height from floor: max $\geq$ 800mm, min $\leq$ 500mm  
 Cushion incline range:  $\geq$ 10°  
 Back-rest incline range: 105° ~170°  
 Chair revolve range: 360°  
 Electric dental chair also available



**SW6306 HIGH SPEED TURBINE HAND PIECE**



**SW6307 DENTAL DEVICE**



**SW6301 DENTAL BRUSH**



**SW6302  
DISPOSABLE DENTAL DEVICE**  
 Saliva ejector, probe, mirror, etc.

# STERILIZER



**YXQ.SG41.280 PORTABLE STEAM PRESSURE STERILIZER**



**YXQ.DY.250 TABLE STEAM PRESSURE STERILIZER (PROGRAM CONTROLLED AND DOOR SELF-LOCKED)**

## LS-B50L VERTICAL PRESSURE STEAM STERILIZER

| ITEM                | YXQ.SG41.280 | YXQ.DY.250    | LS-B50L     |
|---------------------|--------------|---------------|-------------|
| Working Pressure    | 0.15MPa      | 0.22MPa       | 0.14MPa     |
| Working Temperature | 126° C       | 135° C ± 1° C | 126° C      |
| Valid volume        | 16 L         | 20 L          | 50L         |
| Power Supply        | 2KW          | 1.5KW         | 4.5KW       |
| Net Weight          | 14kg         | 36kg          | 82kg        |
| Casing              | 46x44x65cm   | 84x61x49cm    | 55x75x126cm |



## **SJ-C-II(PLASTIC POWDER COATED, FOUR TUBE) DIRECTIONS ON THE ULTRAVIOLET LIGHT STERILIZER**

Power voltage: -220V 50Hz  
Capacity: 60W(2 tubes), 120W(4 tubes)  
Ultraviolet light wave: 253.7nm



## **TS-5L/H, TS-10L/H, TS-20L/H ELECTRIC STAINLESS STEEL DISTILLED WATER EQUIPMENT**

Water output: 5 liter/hour, 10 liter/hour, 20 liter/hour  
Heater power: 4.5kw, 7.5kw, 13.5kw





### SW/KCW-5T ULTRASONIC NEBULIZER

Max capacity of medicine cup: 50ml  
Capacity of water trough: 250ml  
Nebulizing amount: 3ml/min~5ml/min



### SW/KCW-5TYD ULTRASONIC NEBULIZER

Max capacity of medicine cup: 50ml  
Capacity of water trough: 250ml  
Nebulizing amount: 3ml/min~5ml/min  
Adjustable electronic timing function (59 minutes) and parameter LCD display



### BD-300B DISPOSABLE SYRINGE DESTROYER

It can not only melt the needle but also destroy the syringe body conveniently and quickly, only within 5 seconds and protect patients from the used syringe's infection.  
Power: 220V/110V AC or DC  
size: 328x248x205mm



### BF-500+ FETAL DOPPLER

It can detect the fetal heart rate (FHR) as early as 12 weeks of pregnancy to the baby's first moment of life and can be used in hospital, clinic, home etc.  
Power: 9.6V rechargeable batteries  
Ultrasonic frequency: 2.5 MHz  $\pm$  0.25 MHz  
Ultrasonic intensity:  $\leq$  5mW/cm<sup>2</sup>  
Work time: >6H



### SW/YP930 INFANT INCUBATOR

- Air-temp or skin-temp servo-controlled by Micro-computer, set-temp, real-temp and aux-temp display in three windows, with RS232 connector
- Auto-adjust heating power, display heating power 0~100%
- Electronical stands
- With temperature deviation, over-temp, sensor failure, fan low-speed and power broken alarm and so on
- $\geq$ 37°C temperature set, inclination of the bassinet adjustable steeplessly
- Isolated system alarm
- Equipped with: guardrail, whole trough(AL), quiet operation windows, double walls hood, airtemp sensor
- When opening the front door, automatic air curtain come

# SUCTION UNIT



## YX920S ELECTRIC SUCTION APPARATUS

Max. Vacuum:  $\geq 0.09$ MPa (680mmHg)  
 Adjustable vacuum range: 0.02MPa - Max. vacuum  
 Flow rate:  $\geq 20$  L/min  
 Noise:  $\leq 60$ dB(A)  
 Storage bottle: 1000ml



## YB-SXT-1A PORTABLE ABSORB PHLEGM UNIT

Max. vacuum:  $\geq 0.08$ MPa (600mmHg)  
 Adjustable vacuum: 0.02Mpa (150mmHg)-Max. vacuum  
 Flow rate:  $\geq 22$ L/min  
 Noise:  $\leq 60$ dB(A)  
 Storage bottle: 1000ml



## YX932S PORTABLE SUCTION APPARATUS

Max. vacuum:  $\geq 0.09$ MPa (680 mmHg)  
 Adjustable vacuum range: 0.02 MPa~0.09 MPa  
 Flow rate:  $\geq 32$  L/min  
 Jar: 2000ml x1 (standard configuration)



## YX930D ELECTRIC SUCTION APPARATUS

Max Vacuum:  $\geq 0.09$ Mpa (680mmHg)  
 Free air flow:  $\geq 30$ L/min  
 Liquid storage bottle: 2500ml x2  
 Vacuum adjustable range: 0.02~0.09MPa (150~680mmHg)



## YX930D-1A Electric Suction Apparatus

Max. vacuum:  $\geq 0.09$ MPa (680mmHg)  
 Adjustable vacuum Range: 0.02MPa - Max. Vacuum  
 Flow rate:  $\geq 30$ L/Min  
 Storage bottle: 2000ml+1000ml



## YX932D ELECTRIC SUCTION APPARATUS

Max. vacuum:  $\geq 0.09$ MPa (680mmHg)  
 Adjustable vacuum Range: 0.02MPa - Max. Vacuum  
 Flow rate:  $\geq 32$ L/Min  
 Storage bottle: 2000ml x2



## YX980D ELECTRIC SUCTION APPARATUS

Max. vacuum:  $\geq 0.09$ MPa (680 mmHg)  
 Adjustable vacuum range: 0.02 MPa~0.09 MPa  
 Flow rate:  $\geq 80$  L/min  
 Storage bottle: 4000 ml x2



## DY-1A LOW-PRESSURE SUCTION APPARATUS

Max. vacuum:  $\geq 18\sim 20$ kpa  
 Adjustable vacuum: 2kPa~Max. vacuum  
 Flow rate:  $\geq 5$ L/min  
 Storage bottle: 1000ml



# LABORATORY



## SWXSZ-107 207 SERIES BIOLOGICAL MICROSCOPE

| OUTFIT       | MODEL   |                 |                     |                |                 |                     |
|--------------|---|-----------------|---------------------|----------------|-----------------|---------------------|
|              | SW/XSZ<br>-107BN  | SW/XSZ<br>-N107 | SW/XSZ<br>-N107C CD | SW/XSZ<br>-207 | SW/XSZ<br>-N207 | SW/XSZ<br>-N207C CD |
| Viewing Head | Monocular Head, at 45° Rotatable                        |                 |                     |                |                 |                     |
|              | Sliding Binocular Head, at 45° Rotatable                | ●               | ●                   | ●              |                 |                     |
|              | Compensation Free Binocular Head, at 30° Rotatable      |                 |                     |                | ●               | ●                   |
|              | Trinocular Head for Photography                         |                 |                     |                |                 | ●                   |
|              | Trinocular Head for Video Camera (C Mount)              |                 |                     | ●              |                 |                     |
| Objective    | Achromatic DIN 4X, 10X, 40X(S), 100X(S, Oil)            | ●               | ●                   | ●              | ●               | ●                   |
| Eyepiece     | Wide Field Eyepiece WF10X, P16X                         | ●               | ●                   | ●              | ●               | ●                   |
| Stage        | Double Layers Mechanical Stage 140X 140mm/75X75mm       |                 | ●                   | ●              | ●               | ●                   |
|              | Single Layers Mechanical Stage 140X 140mm/75X75mm       | ●               |                     |                |                 |                     |
| Condenser    | Abbe NA1.25 with Iris Diaphragm & Filter                | ●               | ●                   | ●              | ●               | ●                   |
| Illumination | Built-in Illumination, Halogen Lamp 110V or 220V/6V/20W | ●               | ●                   | ●              | ●               | ●                   |
|              | Plan-concave Mirror                                     | ●               | ●                   | ●              | ●               | ●                   |



## SWXSZ-20 SERIES BIOLOGICAL MICROSCOPE SWXSZ-2001/2002

| OUTFIT                          | MODEL  |       |      |       |
|---------------------------------|--|-------|------|-------|
|                                 | 2001   | 2001A | 2002 | 2002A |
| Achromatic Objectives           | 4X, 10X, 40X (S), 100X (S, Oil)  |       |      |       |
| Wide Field Eyepiece             | WF10X, WF16X   |       |      |       |
| Sliding Binocular Head          | Inclined at 45°, 360° Rotatable<br>Interpupillary Distance 55~75mm                 |       |      |       |
| Nosepiece                       | Quadruple  | ●     |      | ●     |
|                                 | Quintuple  |       | ●    | ●     |
| Double layer Mechanical stage   | Size 140x140mm, Moving Range 75x45mm   |       |      |       |
| Focusing                        | Coaxial Coarse and Fine Adjustment, Focusing Range 30mm, Focusing Interval 0.002mm |       |      |       |
| Illumination, condenser NA=1.25 | 6V/20W Halogen Lamp, Adjustable Brightness   | ●     | ●    |       |
|                                 | Kohler Illumination  |       |      | ●     |
| Power                           | AC Input 220V or 110V, DC Output 6V  |       |      |       |



## MONOCULAR BIOLOGICAL MICROSCOPE

| OUTFIT                | MODEL                                  |                |                |                |              |
|-----------------------|--|----------------|----------------|----------------|--------------|
|                       | SW/XSZ<br>-2009                        | SW/XSP<br>-101 | SW/XSP<br>-102 | SW/XSP<br>-103 | SW/SME<br>-F |
| Achromatic Objectives | 4X, 10X, 40X (S)                       | ●              | ●              | ●              | ●            |
|                       | 100X (S, Oil)                          | ●              |                |                |              |
| Eyepiece              | P10X                                   | ●              | ●              | ●              |              |
|                       | H5X, H10X, H16X                        |                | ●              |                |              |
|                       | WF 10x                                 |                |                |                | ●            |
| Nosepiece             | Triple                                 |                | ●              | ●              | ●            |
|                       | Quadruple                              | ●              |                |                |              |
| Stage                 | Plain Stage with Slide Clips           |                | ●              | ●              |              |
|                       | Double Layer Mechanical Stage          | ●              |                |                | ●            |
| Focusing              | Separate Coarse                        |                | ●              | ●              |              |
|                       | Coaxial Coarse                         | ●              | ●              |                | ●            |
| Condenser             | Single Lens with Disc Diaphragm        |                | ●              |                |              |
|                       | Abbe NA1.25 Condenser                  | ●              |                | ●              | ●            |
| Illumination          | Plan-concave Mirror                    |                | ●              |                |              |
|                       | Incandescent Lamp 220V/20W or 110V/20W |                |                | ●              | ●            |
|                       | Halogen Lamp 6V/20W                    | ●              |                |                |              |
| Power                 | AC Input 220V or 110V                  |                |                |                |              |



### SW7101, 7102, 7105 SLIDE

Ground edges  
Ungrounded edges  
Frosted one end, ground edges




### SW7102 COVER SLIP

Size: 18 x 18mm, 20 x 20mm, 22 x 22mm, 24 x 24mm,  
24 x 32mm, 24 x 40mm, 24 x 50mm, 24 x 60mm, etc.



### SW715 URIC TEST STRIPS SW716 BLOOD TEST STRIPS



### SW752 RAPID DIAGNOSTIC TESTS


HCG Midstream, HBsAg, HCV, HIV, Morphine, etc.



### BLOOD LANCET SERIES SW7351 FLAKE TYPE SW7352 TWIST TYPE SW7353 CAPPED TYPE



### SW7355 SAFETY TYPE SW7356 AUTO-SAFETY TYPE



### SW738 VACUUM TUBE

Glass tube or plastic tube available



### SW802 ULTRASONIC GEL

250ml, 2.5L, 5L etc.


# LABORATORY PRODUCTS




**SW730 APPLICATOR**  
Length: 3", 6", 8", etc  
Material: wooden, plastic  
With or without cotton tip




**SW731 PIPETTOR**  
Fixed volume  
Adjustable volume



**SL3 SERIES CENTRIFUGE TUBE**



**SL4 SERIES PIPETTE TIPS**



**SL5 SERIES CULTURE PLATES & MICROPLATES**



**SL6 SERIES PASTEUR PIPETTES**



**SW732 TABLE TYPE CENTRIFUGE**

|            | 80-1      | 80-2      | 90-1      |
|------------|-----------|-----------|-----------|
| MAX. SPEED | 4000r/min | 4000r/min | 5000r/min |
| MAX. RCF.  | 1790Xg    | 2130Xg    | 3130Xg    |
| POWER      | 25w       | 40w       | 40w       |



**HIGH SPEED MICROHEMATOCRIT CENTRIFUGE**

|              | SH120-1                   | SH120-2                 |
|--------------|---------------------------|-------------------------|
| MAX.SPEED    | 12000r/min                | 12000r/min              |
| MAX.RCF.     | 15300Xg                   | 15300Xg                 |
| POWER OUTPUT | 100W                      | 350W                    |
| CAPACITY     | 24PCS CAPILLARY<br>Ø1.5mm | 24PCS CAPILLARY<br>Ø3mm |
| TIME RANGE   | 0~10min                   | 0~99min                 |
| N.WT.        | 11.5KGS                   | 11.5KGS                 |
| REGULATION   | STEP                      | STEPLESS                |



**SW820 MEDICAL PACKING**

- Self-sealing flat pouch
- Heat-sealing flat pouch
- Heat-sealing gusseted pouch
- Heat-sealing flat reel
- Heat-sealing gusseted reel
- Self-sealing connecting pouches
- Heat-sealing connecting pouches



**SW804 DISPOSABLE SURGICAL RAZOR SINGLE OR DOUBLE LAYER**



**SW805 DIAGNOSTIC PENLIGHT**



**SW808 WOODEN TONGUE DEPRESSOR**

- Size: Adult, Child
- Material: wooden, plastic
- Individual sterile packing available



**SW809 SURGICAL BRUSH**



**SW810 DISPOSABLE SURGICAL BRUSH**



**SW812 TEST TUBE BRUSH**



**SW815 DISPOSABLE FORCEPS**